
Insufflator 50L

Instructions for Use **en**

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Insufflator for laparoscopy Type FM134

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en

This instructions for use contains information that is subject to copyright. All rights reserved. This instructions for use should not be photocopied, duplicated on microfilm or otherwise copied or distributed, completely or in part, without the approval of W.O.M. WORLD OF MEDICINE GmbH.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in this instructions for use bear registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or unclarities in this instructions for use.

Copyright © W.O.M. WORLD OF MEDICINE GmbH

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E-mail: info.berlin@novanta.com

CE 0197

CE marking according to Directive 93/42/EEC

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Symbols and description

	Follow instructions for use (white image on a blue background)
	Consult instructions for use
	Medical Device
	Authorized for Sale or use by Physician only
	Caution
	Type BF applied part
	Equipotentiality
	Degrees of protection provided by enclosures (IP-Code)
	Alternating current
	Service
	Waste management
	Catalogue number
	Batch code
	Serial number

	Date of manufacture (YYYY-MM-DD)
	Manufacturer
	Non sterile
	Do not resterilize
	Do not reuse
	Sterilized using ethylene oxide
	Use by date (YYYY-MM-DD)
	Quantity
	Number of autoclaving cycles
	Do not use if package is damaged
	Keep dry
	Top-Bottom
	Fragile
	Keep away from sunlight



	Protect from heat and radioactive sources
	Humidity limitation
	Atmospheric pressure limitation
	Temperature limit
	Storage conditions
	Transport conditions
	ON/OFF (push-push)
	Interface data transfer
	Only to be used with medical-grade CO2
	Complies to IEC 60601-1 including US and Canadian deviations in their respective valid version

Table of Contents

1	Important User Notes	9
2	Safety Information	10
3	General Information	11
3.1	Device Description	11
3.2	Intended Use and Contraindications	11
3.2.1	Intended Use	11
3.2.2	Contraindications	11
3.3	General Warnings and Precautions	12
3.3.1	General Warnings	12
3.3.2	Precautions	15
3.4	Indication-Specific Warnings	17
3.4.1	General Warnings Laparoscopy	17
3.4.2	General Warnings Pediatric	18
3.4.3	General Warnings Bariatric	19
3.4.4	General Warnings Vessel Harvest	22
4	Initial Use of the Device	24
4.1	Scope of Delivery	24
4.2	Preparing the Device	24
4.3	Gas Connection	25
4.3.1	Connecting a Gas Bottle	27
4.3.2	Connecting to Central Gas Supply	28
4.3.3	Switching Device Off	28
5	Operating the Device - General	29
5.1	Front of the Device	29
5.2	Rear of the Device	29
5.3	Monitor with Touch Screen	30
5.4	Switching Device On	31
5.4.1	Insufflation Tube Sets	31
5.4.2	Insufflation Tube Connection	32
5.4.3	Gas Heating	32
5.4.4	Displaying/Selecting Insufflation Operating Mode	33
5.4.5	Setting the Nominal Pressure - General	33
5.4.6	Setting the Nominal Flow - General	34
5.5	Gas Consumption Display	35
5.5.1	Starting/Stopping Insufflation	35
5.5.2	Switching Device Off	35
6	Device Use and Control in Standard Operating Mode	36
6.1	Device-Specific Risks Associated with Use of the Device in Standard Operating Mode	36
6.2	Selecting Standard Operating Mode	37
6.3	Presetting Nominal Pressure in Standard Operating Mode	37
6.4	Presetting Nominal Flow in Standard Operating Mode	38
6.5	Performing the Function Test (in Standard Operating Mode)	39
6.5.1	Filling Tube System with CO ₂	40
6.6	Using the Device during Surgery	40
6.6.1	Insufflating with the Veress Cannula	41
6.6.2	Insufflating with the Trocar	41
6.6.3	Stopping Insufflation	41
7	Device Use and Control with Pediatric Operating Mode	43
7.1	Device-Specific Risks Associated with Use of the Device in Pediatric Operating Mode	43
7.2	Selecting Pediatric Operating Mode	44
7.3	Presetting Nominal Pressure in Pediatric Operating Mode	44
7.4	Presetting Nominal Flow in Pediatric Operating Mode	45
7.5	Performing the Function Test in Pediatric Operating Mode before Using the Device during Surgery	46
7.5.1	Filling Tube System with CO ₂	48
7.6	Using the Device during Surgery	48
7.6.1	Insufflating with the Veress Cannula	48
7.6.2	Insufflating with the Trocar	49
7.6.3	Stopping Insufflation	49
8	Device Use and Control with Bariatric Operating Mode	51

8.1	Device-Specific Risks Associated With Use of the Device in Bariatric Operating Mode	51
8.2	Selecting Bariatric Operating Mode.....	53
8.3	Presetting Nominal Pressure in Bariatric Operating Mode.....	53
8.4	Presetting Nominal Flow in Bariatric Operating Mode	54
8.5	Performing the Function Test in Bariatric Operating Mode before Using the Device during Surgery	54
8.5.1	Filling Tube System with CO ₂	56
8.6	Using the Device during Surgery	56
8.6.1	Insufflating with the Veress Cannula	56
8.6.2	Insufflating with the Trocar.....	57
8.6.3	Stopping Insufflation	57
9	Device Use and Control in Vessel Harvest Operating Mode	59
9.1	Device-Specific Risks Associated with Use of the Device in Vessel Harvest Operating Mode	59
9.2	Selecting Vessel Harvest Operating Mode.....	60
9.3	Presetting Nominal Pressure in Vessel Harvest Operating Mode.....	60
9.4	Presetting Nominal Flow in Vessel Harvest Operating Mode	61
9.5	Performing the Function Test in Vessel Harvest Operating Mode before Using the Device during Surgery	61
9.5.1	Filling Tube System with CO ₂	63
9.6	Using the Device during Surgery	63
9.6.1	Insufflation with Endoscopic Vessel Harvest Instrument	63
9.6.2	Stopping Insufflation	64
10	Configuration Menu (Overview)	65
10.1	Configuration Menu I.....	66
10.1.1	Setting First Nominal Pressure.....	67
10.1.2	Setting the Venting Valve.....	67
10.1.3	Setting the Gas Supply Type	68
10.1.4	Setting the Acoustic Signal Volume	68
10.2	Configuration Menu II.....	69
10.2.1	Setting the Gas Flow Rate.....	69
10.2.2	Setting the Maximum Nominal Pressure.....	69
10.2.3	Setting the Occlusion Warning Signal	70
10.3	Settings Menu	70
10.3.1	Checking Software Version	70
10.3.2	Brightness.....	70
10.3.3	Select Language.....	71
10.3.4	Service Menu	71
11	Safety Functions.....	72
12	Care and Maintenance.....	74
12.1	Cleaning the Hose for the Central Gas Supply.....	74
12.2	Cleaning the Device	74
12.3	Maintenance Intervals	74
12.4	Maintenance by Authorized Service Technician	74
12.5	Replacing the Fuse	75
12.6	Care and Maintenance of Universal Gas Connectors	76
12.6.1	Disassembly	77
12.6.2	Accessories	77
12.7	Care and Maintenance of the High-Pressure Hose	77
12.7.1	Visual Check.....	77
12.7.2	Functional Test.....	77
12.8	Care and Maintenance of Reusable Tube Sets (NOT FOR SALE IN USA)	77
12.9	Cleaning and Disinfecting.....	78
12.10	Sterilization of Reusable Tube Sets (NOT FOR SALE IN USA).....	80
13	Annual Inspection.....	82
13.1	Electrical Safety Test.....	82
13.2	Basic Function Test (in Standard Operating Mode).....	82
13.3	Testing the Pressure Sensors in Standard Mode Operating	83
13.4	Pressure Monitoring Test (in Standard Operation Mode).....	84
13.5	Venting Valve Test	84
13.6	Maximal Device Pressure Test.....	85
13.7	Gas Flow Rate Test	85
14	Electromagnetic compatibility	86
14.1	Electrical Connections.....	86
14.2	Guidance and Manufacturer's Declaration - Electromagnetic Emissions	86

14.3 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity	87
15 Error and Warning Messages	89
16 Technical Data	91
17 Accessories	92
18 Appendix	93
18.1 Test log	93
Glossary	95
Index.....	96

1 Important User Notes

Read the instructions for use carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing or service personnel, or
- damages or malfunction of device and/or accessories.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of its products.

Paragraphs marked with the words WARNING, CAUTION, and NOTE carry special meanings. Sections marked with these words must be given special attention.

Subject to technical changes

Please note

WARNING!

The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third party.

CAUTION!

These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.

NOTE!

These paragraphs contain information to clarify the instructions or provide additional useful information.

2 Safety Information

Federal Law (only for U.S. market)

CAUTION! Federal law restricts this device to sale by or on the order of a physician.

Exclusion of liability

The manufacturer is not liable for indirect, incidental and consequential damages, including, but not limited to, loss of profit. Any liability and applicable warranty becomes null and void if:

- the device and/or the accessories/ peripherals are improperly used, transported, stored, prepared, or maintained;
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories/peripherals;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedules are not adhered to.

The handing over of technical documents does not constitute authorization to make repairs, adjustments or alterations to the device or accessories/ peripherals.

WARNING! Modification of the device Insufflator 50L is not permitted.

Authorized service technician

Only an authorized service technician may perform repairs, adjustments, or alterations on the device or accessories/ peripherals and use the service menu. Any violation will void any applicable warranty. Authorized service technicians are only trained and certified by the manufacturer.

Care and maintenance

The service and maintenance of the device and its accessories/ peripherals has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use. Maintenance of the device may not be performed during surgery or operation of the device.

Contamination


NOTE! Service or maintenance work may not be carried out during surgery.

Before shipping, decontaminate the device and accessories/ peripherals in order to protect the service personnel. Follow the instructions listed in these instructions for use. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to refuse the repair of contaminated devices or accessories/peripherals.

Waste management

The product must be prepared by the operator for disposal. When disposing or recycling the product or its components, make sure to comply with national rules and regulations at all times! A product labeled with this symbol must be disposed separately according to local rules about the waste collection of electrical and electronic devices. The disposal is carried out free by the manufacturer if within the European Union. Please contact your local WOM representative or distributor for additional information about the disposal of your product.

The disposal must be carried out in accordance to RDC 306/2004 or regulations replacing it.

The device Insufflator 50L described in these Instructions for Use corresponds to the FM134 Insufflator.

3 General Information

3.1 Device Description

The device **Insufflator 50L** is a CO₂ insufflator intended to be used during diagnostic and/or therapeutic laparoscopic procedures to create a cavity through insufflation.

The device measures the pressure within the cavities and compares the nominal with the actual pressure in the cavity. The function of this device is to maintain the nominal pressure. Any overpressure within the cavity is lowered to the preset nominal pressure by the automatic venting system.

Essential Performance

The **Insufflator 50L** includes the following applications:

Standard operating mode → Insufflation for adults.

Pediatric operating mode → Insufflation for infants and children (neonatal).

Bariatric operating mode → Insufflation for morbidly obese patients.

Vessel Harvest operating mode → Insufflation for endoscopic vessel harvesting procedure.

- **Standard** operating mode, **Pediatric** operating mode, and **Bariatric** operating mode of the device are used in conjunction with a laparoscope to fill and distend a peritoneal cavity with gas.
- **Vessel Harvest** operating mode is used to create a cavity along the saphenous vein and/or the radial artery during an endoscopic vessel harvesting procedure.

3.2 Intended Use and Contraindications

3.2.1 Intended Use

The device **Insufflator 50L** is a CO₂ insufflator intended for use during diagnostic and/or therapeutic endoscopic procedures to distend a cavity by filling it with gas. The **Standard**, **Pediatric** and **Bariatric** operating modes of the device are indicated to fill and distend a peritoneal cavity with gas during a laparoscopic procedure. The **Pediatric** operating mode is specifically indicated for pediatric laparoscopic procedures. **Vessel Harvest** operating mode is indicated for use during endoscopic vessel harvesting procedures to create a cavity along the saphenous vein or radial artery.

3.2.2 Contraindications

The device may not be used to fill an abdomen with CO₂ if a laparoscopy is contraindicated. Please consult the manual of your laparoscope for absolute and relative contraindications. The device is not suitable for hysteroscopic insufflation procedures, i.e., it may not be used to distend the uterus.

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg.

The device may not be used for the endoscopic harvesting of vessels if this surgical application is contraindicated. Please consult the manual of the instrument for absolute and relative contraindications.

3.3 General Warnings and Precautions

3.3.1 General Warnings

WARNING!

Condensation / Water penetration

Protect device from moisture. Do not use if moisture has penetrated the device.

WARNING!

Original accessories

For your own safety and that of your patient, use only original accessories.

WARNING!

Check all factory settings

Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

WARNING!

Technique and procedures

Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

WARNING!

Not explosion-proof

The device is not explosion-proof. Do not operate the device in the vicinity of explosive anesthetic gases and not in the vicinity of oxygen-enriched environments.

WARNING!

Risk of electrical shock

To prevent electrical shock, do not open this device. Never open this device yourself. Notify the authorized service technicians of any required repairs.

WARNING!

Risk of electrical shock

To avoid the risk of electrical shock, only use this device when connected to a properly grounded power supply network.

WARNING!

Replace the fuse

Replace the fuse only with a fuse of the same type and rating (see Chapter Accessories [► 92]).

WARNING!**Professional qualification**

The instructions for use do not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualifications working under the direction and supervision of a physician.

**WARNING!****Functional test**

The functional test must be performed by the user prior to each surgery.

**WARNING!****Sterile media and accessories**

Always work exclusively with sterile substances and media, sterile fluids, and sterile accessories if so indicated.

**WARNING!****Reprocessing of sterile disposable products**

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

**WARNING!****Replacement device and accessories**

In case the device or any of the accessories fail during surgery, a replacement device and replacement accessories should be kept within close proximity to be able to finish the operation with the replacement components.

**WARNING!****Device defect**

If a device defect is suspected or confirmed, do not use it. Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!****Positioning the patient**

Always position the patient lower than the device to prevent body fluids from leaking into the insufflation tube. Actual pressure may increase and fluid may penetrate the insufflation tube if the patient is repositioned during surgery. If this occurs, immediately disconnect the insufflation tube. When the patient is repositioned onto his or her side, internal tissue may block the insufflation channel. Always insufflate through the elevated side of the patient.

**WARNING!****Removing the insufflation tube**

Always disconnect the insufflation tube after ending surgery and before switching off the device to prevent backflow of bodily fluids. Fluid may penetrate the insufflation tube whenever you change the gas bottle and/or when you stop the gas flow during the operation. If this happens, you must immediately disconnect the insufflation tube from the trocar or from the device.

**WARNING!****Backflow**

Body secretions or contaminated gas may backflow into the device through the insufflation tube if

- a filter is not used,
- the actual pressure is higher than the nominal pressure or
- the automatic venting valve is activated.

**WARNING!****Gas flow**

A high gas flow can occur due to large leaks within the surgical system or instrument. This can result in a false actual pressure reading, which in turn may endanger the patient. In case of a disrupted gas flow, you should therefore inspect device, tube, and instruments immediately. Surgical procedures should be performed with a gas flow of 4 to 10 l/min. An even lower gas flow is recommended for diagnostic purposes.

**WARNING!****Keep filled CO₂ bottle on hand**

Always keep a filled CO₂ bottle on hand ready for replacement. This avoids having to interrupt surgery due to a lack of insufflation gas (see Chapter Gas Connection [► 25]).

**WARNING!****Contamination**

Do not use device and/or accessories if signs of contamination are detected. Make sure the device or/and accessories can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!****Fatigue symptoms**

When there is a high level of CO₂ consumption, you should make sure to supply the operating area with enough fresh air, since an increasing CO₂ level in the air can cause the medical personnel to suffer fatigue symptoms, an inability to concentrate, unconsciousness, or even death.

**WARNING!****Limited venting system**

The venting rate of the automatic venting system is limited. Always monitor the actual pressure when using additional insufflation sources.

WARNING!**Contaminated filter**

Replace a contaminated filter immediately during surgery to ensure unhindered gas flow.

**WARNING!****Connecting the tube**

Always use the proper tube set for the device. The tube outlet may only be connected to instruments which are intended for intra-abdominal CO₂ insufflation.

**WARNING!****Electronic device control**

Do not close the valve (stopcock) at the trocar shaft during surgery. The electronic control unit of the device adjusts the actual pressure as desired.

**WARNING!****Medically pure CO₂**

Make sure to use only medically pure CO₂. Other gases (e.g., helium, N₂O, argon), mixtures of gases, gases with entrapped liquids, or polluted gases must not be used with this device.

**WARNING!****Increased airway pressure/compression of the vena cava**

When using the device on children, an increased risk of elevated airway pressure and/or compression of the vena cava (low input syndrome) exists.

**WARNING!****Emphysema**

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries (> 200 min.), the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar accesses immediately.

**3.3.2 Precautions****CAUTION!****In case of serious incident**

Please report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the State in which the user and/or patient is established.

**CAUTION!****Cleaning the Device**

Do not sterilize the device.

**CAUTION!****Service connection**

Access to the service menu is restricted to authorized service personnel. The connected equipment must comply with the standard EN 60950 in the currently valid version. Do not connect a device to the service connection during surgery.

**CAUTION!****Peripheral Devices**

Additional peripheral equipment connected to interfaces of the medical monitor has to meet the requirements of the following specifications in the respective current valid version: IEC 60601-2-18 / EN 60601-2-18 for endoscopic devices and IEC 60601-1 / EN 60601-1 for electrical medical devices. All configurations have to comply with IEC 60601-1 / EN 60601-1 specifications. Whoever connects additional equipment to signal output or signal input is considered the system configurator and as such is responsible for complying with requirements of the standard IEC 60601-1 / EN 60601-1.

**CAUTION!****Incorrect voltage**

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

**CAUTION!****Endoscope**

The device may only be connected with endoscopes designed for and featuring the technical specification permitting such a combined use. Any utilized endoscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600. Combining/connecting with other devices generates a medical electrical system (MES). The system configurator is responsible for compliance with the standard IEC 60601-1 / EN 60601-1 in its latest version.

**CAUTION!****Electrical interference**

(See Chapter Electromagnetic compatibility [▶ 86]). Care was taken during the development and testing of this device that electrical interference of or from other devices or instruments was practically eliminated. However, if you still detect or suspect such interference, please follow these suggestions:

- Move this, the other, or both devices to a different location
- Increase distance between used devices
- Consult an electro-medical expert

**CAUTION!****Use of other accessories, other transducers and cables**

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM (see Chapter Glossary) as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

CAUTION!**ME Device in Rack**

The ME device may not be used when stacked or stored directly adjacent to or with other devices. If it should be necessary to operate the device close to stacked with other devices, the ME device or ME system should be monitored to ensure it works properly as configured.

**3.4 Indication-Specific Warnings****3.4.1 General Warnings Laparoscopy****WARNING!****Idiosyncratic reactions**

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

**WARNING!****CO₂ absorption**

CO₂ is absorbed during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

**WARNING!****Metabolic and cardiac reactions**

Insufflating CO₂ may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solutions that are not preheated to body temperature.

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO₂ embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**WARNING!****Additional insufflation sources / automatic venting system**

Make sure the automatic venting system is activated (see Chapter Configuration Menu (Overview) [▶ 65]). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

3.4.2 General Warnings Pediatric

**WARNING!****Hypercapnia**

Because pediatric patients are especially susceptible to hypercapnia, it is recommended to establish an end-tidal CO₂ monitoring routine.

**WARNING!****Procedures with children**

Only those who are specially trained and qualified for procedures with children or endoscopic vessel harvesting procedures may use this device for these purposes.

**WARNING!****Gas flow limit**

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg (approximately 55 US pounds).

**WARNING!****Recommended work settings**

The flow values listed for laparoscopic procedures performed on newborns, infants, and children are only suggested values. The selection of the suitable flow and pressure values is solely the responsibility of the attending physician.

**WARNING!****Pneumolabium/pneumoscrotum**

Children are at risk of a pneumolabium or pneumoscrotum.

WARNING!**Increased airway pressure**

When laparoscopic procedures are performed on children, the increased intra-abdominal pressure also increases the risk for higher airway pressures. Always strictly monitor respiration and airway function when performing laparoscopic procedures on children younger than 12 years of age.

**WARNING!****Compression of the vena cava**

When insufflating the abdomen of a child with medical CO₂, an increased risk of compressing the vena cava exists. This risk can be reduced by monitoring the systolic and diastolic blood pressure during the entire surgery.

**WARNING!****Haemodynamic stability**

A laparoscopy performed on children younger than 12 years of age may result in problems of the haemodynamic system due to the CO₂ content in the blood. It is recommended to increase the breathing rate of the patient and to work with low flow values and pressure values not exceeding 12 mmHg. The patient's circulatory system should be monitored at all times.

**WARNING!****Hypothermia**

The insufflation gas flow usually drops significantly after the target pressure has been reached and it is then only required to maintain the abdominal pressure. However, leaks within the abdomen or the instrument can lead to a constant gas flow of above 1 l/min. When operating on children younger than 12, a gas flow of more than 1 l/min poses an increased risk of hypothermia for the patient. Corresponding measures to prevent hypothermia include the use of blankets or pre-warmed gas. The patient's body temperature has to be monitored at all times during surgery.

3.4.3 General Warnings Bariatric

WARNING!**Altered Respiratory Physiology**

Always monitor the patient's respiratory functions during the entire surgery. The larger body mass supported by the thoracic cage and the larger amount of fat in the abdominal cavity may reduce the elasticity of the thoracic wall. In addition, the increased intra-abdominal pressure secondary to insufflation may alter the normal physiological lung parameters thus resulting in a reduction of the functional lung volume. Shallow, rapid breathing is symptomatic of this condition. Even modest physical stress causes a tremendous increased demand for oxygen, which stands in contrast to the ineffective respiratory musculature that requires more oxygen because it must overcome the reduced elasticity of the thoracic cage. The functional capacity of the lungs is small and even moderate stress can lead to respiratory failure.

**WARNING!****Subcutaneous emphysema**

When puncturing the thicker abdominal wall of morbidly obese patients with the Veress cannula or the trocar, carefully monitor the correct position of the instrument in the abdomen.

**WARNING!****Idiosyncratic reactions**

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

**WARNING!****CO₂ absorption**

CO₂ is absorbed during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

**WARNING!****CO₂ supersaturation**

To avoid generating CO₂ supersaturation, an increased level of respiratory activity is required. An overweight patient's oxygen demand and carbon dioxide production are greater and increase faster under physical stress than do those of patients with normal body weight.

**WARNING!****Heart and cardiovascular insufficiency**

Constantly monitor all heart and cardiovascular parameters during surgery since morbidly obese patients have an increased risk of heart and cardiovascular insufficiencies.

**WARNING!****Metabolic and cardiac reactions**

Insufflating CO₂ may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

WARNING!**Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solutions that are not preheated to body temperature.

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO₂ embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**WARNING!****Emphysema**

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries (> 200 min.), the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar accesses immediately.

**WARNING!****Additional insufflation sources / automatic venting system**

Make sure the automatic venting system is activated (see Chapter Configuration Menu (Overview) [▶ 65]). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

3.4.4 General Warnings Vessel Harvest

WARNING!

Procedures with children

Only those who are specially trained and qualified for procedures with children or endoscopic vessel harvesting procedures may use this device for these purposes.

WARNING!

Instrument used for CO₂ insufflation

Before using the insufflator to endoscopically harvest vessels, please check whether the instrument used is intended for CO₂ insufflation.

WARNING!

Pneumoperitoneum

When a vessel is harvested from the leg of a patient with a perforated groin, it is possible for CO₂ to reach the abdomen and cause a pneumoperitoneum. Make sure the abdomen does not fill with CO₂ during surgery.

WARNING!

Idiosyncratic reactions

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

WARNING!

CO₂ absorption

Due to the special surgical procedures - start of the heart bypass operation, and the endoscopic removal of the vessel - special care has to be taken as CO₂ is always absorbed through the tissue of the patient during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption.

WARNING!

Metabolic and cardiac reactions

Due to the special surgical conditions - start of the heart bypass surgery and vessel harvesting - it is especially important to remember the existing risk of metabolic acidosis when insufflating with CO₂. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

WARNING!**Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

WARNING!**Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO₂ embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

Delivery inspection**Returning the device****Setting up****4 Initial Use of the Device**

Always check all parts and accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized service company.

4.1 Scope of Delivery

- Device **Insufflator 50L**
- Instructions for Use
- Power supply cord

Use the original packaging material if the device has to be returned. The manufacturer is not responsible for transport damages due to insufficient or unsuitable packaging.

- Name of owner
- Address of owner
- Device type and model
- Serial number of the equipment (see identification plate)
- Detailed description of defect

4.2 Preparing the Device

Place the device on a flat surface free of vibration located in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in Chapter Technical Data [► 91].

WARNING!**Not explosion-proof**

The device is not explosion-proof. Do not operate the device in the vicinity of explosive anesthetic gases and not in the vicinity of oxygen-enriched environments.

WARNING!**Position of the device**

Position the device in such a way that it is easy to operate and switch off.

CAUTION!**ME systems**

The medical electrical (ME) device is suitable for integration in ME equipment systems (see Chapter Glossary). Operation of the ME device in vicinity of non-ME devices may result in voiding the intended use of the ME device.

CAUTION!**Device setup**

Device should be positioned outside of the sterile area in such a way that

- it can be easily disconnected,
- it is easy to use and switch off and on,
- it allows an easy monitoring of the display values, device functions, and access to the control elements.

CAUTION!**Errors and malfunctions of the device**

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

Make sure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The mains connection cable may be plugged only into a properly installed, grounded safety wall socket (shockproof socket) (see DIN VDE 0100-710). Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

The power connection must be equipped with a grounding contact. Use the original power cable (if included in scope of delivery) to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.


Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 and IEC 60320-C13. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

Integrate the device into the potential equalization system as specified by local safety rules and regulations.

Mains connection**Grounding contact****Only for U.S. operators****Potential equalization****4.3 Gas Connection****WARNING!****Medically pure CO₂**

Make sure to use only medically pure CO₂. Other gases (e.g., helium, N₂O, argon), mixtures of gases, gases with entrapped liquids, or polluted gases must not be used with this device.

Use a high-pressure hose to connect a CO₂ gas bottle to the gas supply connection or use a low pressure hose to connect to a centralized CO₂ gas supply.


The device is equipped with a universal gas connector that can be configured with different adapters to support both central gas and bottle gas supply.

Universal gas connector**Features**

- Angled shape to reduce the depth needed for the installation of the device
- Exchangeable sinter filter for the gas inlet to protect the insufflator from dirt particles
- Maximum gas pressure: 80 bar
- Pre installed US-connector for bottle gas supply

Exchangeable sinter filter


The gas connector has an exchangeable sinter filter to protect the insufflator from dirt particles. This filter has to be checked during the annual test and replaced at least every two years. During the maintenance the sinter filter has to be exchanged. The sinter filter is available as an accessory (see Chapter Accessories [► 92]). The sinter filter may be blocked or soiled due to unclean gas or gas supply residues. Heavy soiling is especially noticeable during higher flow rates. As a result, the actual gas flow is falling short of the nominal gas flow during active insufflation or the device is showing gas supply warnings although gas supply is sufficient. In this case the sinter filter needs to be exchanged.

1. Remove the connector ① using two open-end wrenches size 17 mm and 19 mm.
2. Unscrew the sinter filter ② using a screwdriver if necessary.
3. Screw the new sinter filter finger-tight and fasten the connector to the device.

Angled gas connector

Designation
Universal gas connector CO2 angled - standard

An angled connector is available to reduce the depth needed for the installation of the device. This connector includes the exchangeable sinter filter and allows the use of the different gas supply adapters described below as well.

4.3.1 Connecting a Gas Bottle

CAUTION!

High pressure hose

Always use a high pressure hose to connect gas bottle and device.

The gas bottle must be in a vertical position. The gas bottle pressure may not exceed 80 bar or be less than 15 bar.

CAUTION!

Gas bottles with riser pipe

Gas bottles with riser pipe can release dirt and oily fluids into the device. Do not use a gas bottle with riser pipe.

High Pressure Hose Types

Designation
Hose for CO2 gas bottle DIN:US, 1.5 m
Hose for CO2 gas bottle ISO:US, 1.5 m
Hose for CO2 gas bottle PIN:US, 1.5 m
Hose for CO2 gas bottle CGA:US, 1.5 m

Installation

High Pressure Hose with Connector PIN

1. Connect or remove hose from the device using the open-end wrench size 14 millimeter.
2. Connect or remove hose from the gas bottle with your hand.

High Pressure Hose with Connector DIN

1. Connect or remove hose from the device using the open-end wrench size 14 millimeter.
2. Connect or remove hose from the gas bottle using the open-end wrench size 30 millimeter.

High Pressure Hose with Connector ISO

1. Connect or remove hose from the device using the open-end wrench size 14 millimeter.
2. Connect or remove hose from the gas bottle using the open-end wrench size 32 millimeter.

High Pressure Hose with Connector CGA

1. Connect or remove hose from the device using the open-end wrench size 14 millimeter.
2. Connect or remove hose from the gas bottle using the open-end wrench size 1-1/8 inch.

Switching valve

Designation
Switch valve for insufflators

There is a switching valve available that allows the connection of two gas bottles to the insufflator. This ensures uninterrupted insufflation particularly during the change of the gas bottle.

4.3.2 Connecting to Central Gas Supply

The following low pressure hoses are available for connecting to the central gas supply:

Low pressure hoses

Designation
Hose for CO2 central gas DIN:US, 3 m
Hose for CO2 central gas DIN:US, 5 m
Hose for CO2 central gas NF:US, 3 m
Hose for CO2 central gas NF:US, 5 m
Hose for CO2 central gas UNI:US, 3 m
Hose for CO2 central gas UNI:US, 5 m
Hose for CO2 central gas AGA:US, 3 m
Hose for CO2 central gas AGA:US, 5 m
Hose for CO2 central gas DISS:DISS, 3 m

Installation

- Connect or remove hose from the device using the open-end wrench size 14 mm.
- Connect or remove hose from the central gas supply receptacle on the wall pulling the plug in or out.

An adapter and a low pressure tube is required for connecting to the central gas supply with DISS connection:

DISS adapter for universal gas connector

Hose for CO₂ central gas DISS:DISS

Exchange of the gas supply adapter


Remove the connector ① using two open-end wrenches size 17 mm and 19 mm.

Check if the sinter filter ② needs a change (residues or blocked pores).

If necessary, exchange the sinter filter using a screwdriver.

Install the new connector ① using the open-end wrenches.

Set in the user menu the type of gas supply: central gas or bottle gas supply (see Chapter Configuration Menu (Overview) [► 65]).

4.3.3 Switching Device Off

Disconnect device from power supply

NOTE!

Safe disconnection from power supply


A safe and all-pole disconnection of the product from the power supply is only ensured by disconnecting the power cord.

1. Switch device off: press the ON/OFF button (push-push) ① (Fig. Device front [► 29]).
2. Remove the power cord from the device mains socket ⑨ (Fig. Device rear [► 29]).

The operation of the device is safely terminated.

5 Operating the Device - General

5.1 Front of the Device


Familiarize yourself with the control and function elements at the front of the device.

Fig. 5-1 Device front

- ① ON/OFF button (push-push)
- ② Touch screen display
- ③ Insufflation tube connection
- ④ Gas heating connection for reusable heating tubes
(NOT FOR SALE IN USA)
- ⑤ Gas heating connection for disposable heating tubes

5.2 Rear of the Device

Familiarize yourself with the connection elements at the rear of the device.

Fig. 5-2 Device rear

- ⑥ Service interface
- ⑦ Connection for potential equalization
- ⑧ Fuse holder
- ⑨ Device mains socket
- ⑩ Gas supply connection


Fig. 5-3 Screen display

- ① Insufflation mode selection
- ② Insufflation mode display
- ③ Status display/error and warning messages
- ④ Menu
- ⑤ Increase nominal gas flow
- ⑥ Nominal gas flow display
- ⑦ Decrease nominal gas flow
- ⑧ Actual gas flow display
- ⑨ Start/Stop function field
- ⑩ Gas heater connected/ready
- ⑪ Gas supply display
- ⑫ Gas consumption display
- ⑬ Function field to reset the gas consumption display (re-set button)
- ⑭ Decrease nominal pressure
- ⑮ Nominal pressure display
- ⑯ Increase nominal pressure
- ⑰ Abdominal pressure display

Gas supply displays

Gas supply with gas bottle

5.3 Monitor with Touch Screen

The above depiction of the touch screen also shows all display and function fields.

The status of the gas supply is monitored by the device and indicated with symbols and acoustic signals (see Chapter Safety Functions [▶ 72]) for gas pressure display information).

The following gas bottle pressures are displayed:

	> 50 bar
	40 – 50 bar
	30 – 40 bar
	15 - 30 bar; three acoustics warning signals (beeps) can be heard and the message Change gas bottle is displayed. User is advised to prepare for changing the gas bottle.
	< 15 bar; three warning signals can be heard and the message Check gas supply is displayed. Immediately replace the gas bottle.

If the gas bottle pressure keeps falling, warnings will appear to remind the user to replace the gas bottle immediately. At < 5 bar and again at 0 bar, 5 acoustic signals will sound and the display depicts **Check Gas Supply**. Insufflation stops at 0 bar.

The following house gas supply pressures are displayed:

House gas supply pressure OK

(green)


House gas supply pressure too low

(red)

House gas supply

5.4 Switching Device On

1. Connect the device with the mains power.
2. Connect the gas supply to the gas connection port (see Chapter Gas Connection [► 25]) and open the gas supply.
3. Press the ON/OFF button (push-push). The device turns on.
4. The device performs a device check after being switched on. The touch screen depicts **Device Test** → **Device OK** for 3 seconds after the successful device test. In case the device check fails and an error message is displayed, please see Chapter Error and Warning Messages [► 89] for additional information.
5. The display depicts an overview of the insufflation operating modes. Press the respective function field to select the corresponding desired operating mode (e.g. **Standard**).


5.4.1 Insufflation Tube Sets

Different insufflation tube set types can be connected to the insufflation tube connection at the front of the device.

Art. No.	Description
Disposable Tube Sets	
Z0293	Insufflation tube set with ISO connection, disposable, sterile
Z0260	Heating tube set with ISO connection, disposable, sterile
Reusable Tube Sets (NOT FOR SALE IN USA)	
Z0452-01	Insufflation tube set ISO connector, reusable 20 times, autoclavable
Z0280	Heater tube, ISO connection, reusable 100 times, autoclavable

Disposable tube sets have been sterilized with ethylene oxide according to the validated methods ISO 11135-1 and ISO 10993-7.

5.4.2 Insufflation Tube Connection

1. Always connect a hydrophobic filter to the reusable tube (**NOT FOR SALE IN USA**) or use the disposable tube set with integrated filter.
2. Connect the insufflation tube to the insufflation tube connection and to the Veress cannula/trocar.

ISO 5356-1 connection

A male ISO connector is located at the device (ISO-M). Connect the female ISO connector (ISO-F) of the filter or tube firmly with the ISO-M connector.

5.4.3 Gas Heating

Use the gas heater and the heating tube to insufflate lukewarm gas (37 °C). Connect insufflation tube ① and heater plug ② as depicted in Fig. Connecting the gas heater [► 32].


Fig. 5-4 Connecting the gas heater

- ① Insufflation tube
- ② Gas heating connection for reusable heating tubes (**NOT FOR SALE IN USA**)

CAUTION!

Direct heat

Do not subject the heating tube to direct heat (e.g., endoscope connected with light source) or high room temperatures.

1. Connect the filter and the heating tube.
2. Connect the plug of the heater tube with the appropriate gas heating connection.

NOTE!

Mix-up of gas heater connections

A mix-up cannot occur because of the different shapes of gas heater connections for disposable tube sets and reusable tube sets (**NOT FOR SALE IN USA**).

Once the heating tube has been connected, the status bar ③ (see Fig. Screen display [► 30]) displays **Gas Heater OK** message for 3 seconds.

NOTE!

Prevention of overheating

In order to reduce the risk of overheating, the device switches the heating function off at very low flow rates.

Start gas heating: Press the **Start/Stop** function field. The gas is automatically heated. The gas is automatically heated to 37 °C in approximately 10 minutes.

The message **Gas Heater defective** (3) (see Fig. Screen display [► 30]) displays is displayed for 2 seconds in case of a function defect of the gas heater (e.g., cable flawed or broken, loose plug).

Check gas heating using a different tube. Should the error message be displayed again, you can continue using the device without gas heating by observing the risks for hypothermia. Check gas heating after surgery using a different tube.

If the message **Gas Heater defective → Call service** persists, call an authorized service technician to check the device.

The status bar (3) (see Fig. Screen display [► 30]) displays depicts the following displays in case of excess temperature: **Temperature > 42 °C** (for 3 seconds).

Gas heating ON/OFF

Gas heating error

Gas temperature exceeds 42 °C

WARNING!

Gas temperature > 42 °C

If the temperature sensor measures a gas temperature > 42 °C, gas heating is deactivated. Hot gas in the abdomen can lead to serious injuries.

A warning signal is audible. Gas heating is deactivated.

1. Continue operation, the device insufflates without heating the gas.
2. Check gas heating after surgery using a different tube. Turn device off and on again. Gas heating is reactivated.
3. Should the error message be displayed again, you can continue using the device without gas heating by observing the risks for hypothermia.
4. Call an authorized service technician to check the gas heating.

5.4.4 Displaying/Selecting Insufflation Operating Mode

1. Four device configurations are available:


Standard operating mode → Insufflation for adults.

Pediatric operating mode → Insufflation for infants and children (neonatal).

Bariatric operating mode → Insufflation for morbidly obese patients.


Vessel Harvest operating mode → Insufflation for endoscopic vessel harvesting procedures.

2. Control or selection of the insufflation modes can be done only while insufflation is stopped. Press the function field for 2 s to change the insufflation mode.
3. The display depicts an overview of the insufflation modes. Unavailable modes are depicted in gray and cannot be selected. Press the respective function field to select the corresponding desired mode (e.g. Pediatric).
4. After pressing the function field for the insufflation mode, the display depicts the selected procedure (e.g. Pediatric mode). The displayed parameters correspond with the factory settings or the values set in the configuration menu (see Chapter Configuration Menu I [► 66]).

5.4.5 Setting the Nominal Pressure - General

Press the function fields ▲ or ▼ on the display below the text line **Set Pressure** to set the nominal pressure.


Applies to all modes:

- The nominal pressure can be increased or decreased during insufflation or while insufflation is stopped.
- With every touch of the pressure function field ▲ or ▼ the nominal pressure is increased or reduced in steps of 1 mmHg. When pressing the function field ▲ or ▼ longer than 1.5 seconds, scrolling (in 1 mmHg steps) is activated but only up to the respective, application-dependent pressure safety limit. Once this limit has been reached, the status line depicts the message **Safety limit** and the actual value flashes. To allow settings above the pressure safety limit release the pressure function field for 2 seconds. The display switches back to the status display. Now it is possible to increase the pressure above the pressure safety limit but only to the next pressure safety limit (as far as existent) or to the maximum nominal pressure, preselected in the configuration menu (see Chapter Configuration Menu I [► 66]).
- The latest working settings are stored with interruption of the insufflation, but only if the last value was below the lowest pressure safety limit.
- **Standard** operating mode 15 mmHg
- **Bariatric** operating mode 15 mmHg
- **Pediatric** operating mode between 12 mmHg and 15 mmHg
- **Vessel Harvest** operating mode between 12 mmHg and 15 mmHg
- Otherwise, the set pressure value is reset to the lowest pressure safety limit upon restarting the insufflation.
- The latest settings are reset to the "First nominal pressure" values preselected in the **Configuration menu** when deactivating the device or when changing into the **Configuration menu** or the Mode Selection menu. As soon as the operating mode has been selected, the nominal values are displayed as stored in the configuration menu (see Chapter Configuration Menu I [► 66]).

5.4.6 Setting the Nominal Flow - General

Press the function fields ▲ or ▼ on the display below the text line **Set Flow** to set the nominal flow.

Applies to all modes:

- The nominal flow can be increased or decreased during insufflation or while insufflation is stopped.
- With every touch of the flow function field ▲ or ▼ the nominal flow is increased or reduced in steps of 1 l/min. When working within the range of 0.1 to 2 l/min in **Pediatric** operating mode, the value is increased/reduced 0.1 l/min. When working within the range of 1 to 5 l/min in **Vessel Harvest** operating mode, the value is increased/reduced 0.5 l/min. When pressing the flow function fields ▲ or ▼ longer than 1.5 seconds, fast switching is activated but only up to the next, flow Level, preselected in "gas flow rate" in the Configuration menu (not active for pediatric mode). Changing of nominal flow value by pressing continuously the function field ▲ is also interrupted at flow safety limit, if **Pediatric** or **Vessel Harvest** operating mode is selected and if "flow safety limit" is switched "ON" in Configuration menu.

Once this flow safety limit has been reached, the status line depicts the message **Flow safety limit** and the actual value flashes. To allow settings above the flow safety limit release the function field for 2 seconds. The display switches back to the status display. Now it is possible to increase the flow above the flow safety limit but only to the maximum setting (see definition of fixed maximum gas flow values in Chapter Technical Data [► 91]).

- The last set values are saved after the interruption of insufflation, but only if the last value was below the flow safety limit (5 l/min for **Pediatric** operating mode and 6 l/min for **Vessel Harvest** operating mode). Otherwise the set flow value is reset to the flow safety limit upon deactivation of insufflation.
- The last set values are reset to the Level 1 flow value preselected in "gas flow rate" in the Configuration menu, respectively to the "first nominal gas flow" (see Chapter Configuration Menu (Overview) [► 65]), if **Pediatric** operating mode is selected, when deactivating the device or when changing into the configuration menu or

the **Mode Selection menu**. Upon choice of the respective working mode the nominal settings as stored in the **Configuration menu** are displayed (see Chapter Configuration Menu (Overview) [► 65]).

5.5 Gas Consumption Display

The gas consumption display indicates the insufflated volume of CO₂ in liters since the last resetting of the display. The display indicates values between 0.0 liter and 9.9 liter in increments of 0.1 and between 10 liter to 999 liter in increments of 1.

The gas consumption display can be reset to 0.0 while insufflation is started as well as stopped. Press the gas consumption display/function field to reset the gas consumption display.

Reset the display

5.5.1 Starting/Stopping Insufflation

The **Start** or **Stop** function field is displayed depending on the insufflation status.

Start insufflation:

The **Start** function field is displayed while insufflation is stopped. Press this field to start insufflation.

The status line depicts the selected insufflation mode (e.g. **Mode: Standard**). The insufflation operating mode display depicts **Change Mode**.

Start insufflation by pressing the **Start** function field.

Insufflation started:

The status line depicts **Insufflation** and the type of selected insufflation mode (e.g. **Insufflation: Standard or Veress**).

The "Select Insufflation Mode" and "Menu" keys are depicted in light gray. This signals that these function fields are inactive during this state.

The insufflation operating mode display depicts the selected mode (e.g. **Mode: Standard**).

Stop insufflation:

The **Stop** function field is displayed while insufflation is active.

Stop insufflation by pressing the **Stop** function field.

The status line depicts **Insufflation stopped**.

5.5.2 Switching Device Off

Disconnect device from power supply

NOTE!

Safe disconnection from power supply

A safe and all-pole disconnection of the product from the power supply is only ensured by disconnecting the power cord.

1. Switch device off: press the ON/OFF button (push-push) ① (Fig. Device front [► 29]).
2. Remove the power cord from the device mains socket ⑨ (Fig. Device rear [► 29]).

The operation of the device is safely terminated.

en	
Patients	
Characteristics	
Insufflation parameters	

6 Device Use and Control in Standard Operating Mode

Standard operating mode uses supplied gas to distend the peritoneal cavity during laparoscopic procedures.

Standard operating mode is designed explicitly for laparoscopies performed on normal weight and slightly obese (BMI < 30 kg/m²) patients over the age of 14.

Standard operating mode is suitable for any size of cavity.

The insufflator limits the pressure to maximal 30 mmHg and the gas flow rate to maximal 40 l/min.

The device measures the pressure within the abdomen at short intervals and compares the nominal with the actual abdominal pressure. The function of this device is to maintain the nominal pressure. Any overpressure within the abdomen is lowered to the preset nominal pressure by the automatic venting system.

6.1 Device-Specific Risks Associated with Use of the Device in Standard Operating Mode

WARNING!

Idiosyncratic reactions

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

WARNING!

CO₂ absorption

CO₂ is absorbed during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

WARNING!

Metabolic and cardiac reactions

Insufflating CO₂ may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

WARNING!**Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solution not pre-heated to body temperature.

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO₂ embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.


**WARNING!****Additional insufflation sources / automatic venting system**

Make sure the automatic venting system is activated (see Chapter Configuration Menu (Overview) [► 65]). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.


6.2 Selecting Standard Operating Mode

1. Selection of the insufflation operating mode can be done only while insufflation is stopped. If Standard operating mode is not yet set, press the function field **Insufflation Mode Selection** for 2 seconds to load the insufflation overview.
2. The display depicts the different insufflation operating modes that can be activated. Press the function field **Standard** to select the desired insufflation operating mode.
3. The procedure profile is depicted on the display. The displayed parameters correspond with the factory settings or the values set in the configuration menu (see Chapter Configuration Menu (Overview) [► 65]).

6.3 Presetting Nominal Pressure in Standard Operating Mode

The nominal pressure can be set during insufflation or while insufflation is stopped. Values may range from 1 to maximal 30 mmHg or correspond with the value set in the configuration menu in this area.

Increasing/decreasing nominal pressure:

Briefly press the ▲ or ▼ function field to increase or decrease the pressure. Pressing the ▲ or ▼ function field longer than 1.5 seconds activates scrolling in increments of 1.

Pressure safety limit:

Increasing the nominal pressure displays **Safety limit** in the status line once the value 15 mmHg has been reached. The nominal pressure value 15 mmHg is a safety threshold value. This is where the recommended range for the intra-abdominal pressure ends. Pressing the nominal pressure function field ▲ again does not increase the pressure any further. Release the function field for 2 seconds. Now you can set a value up to 30 mmHg.

CAUTION!

Exceeding pressure safety limit

Exceeding this pressure safety limit is to be decided by and the responsibility of the user/ operator.

6.4 Presetting Nominal Flow in Standard Operating Mode

The nominal flow can be set during insufflation or while insufflation is stopped.

1. Briefly press the ▲ or ▼ function field to increase or decrease the flow. The nominal flow can be increased from 1 to maximal 40 l/min.
 - Briefly press the corresponding field to set values in increments of 1.
 - Pressing the ▲ or ▼ function field longer than 1.5 seconds activates scrolling through the gas flow rates 3, 20, 40 l/min or 40, 20, 3 l/min.
2. The preset values of the configuration menu can be changed individually (see Chapter Configuration Menu I [► 66]). Select a nominal gas flow between 1 and 40 l/min. The preset value is displayed in the nominal flow display. The values for the nominal gas flow refer to a device without connected tube, filter, or instrument. Tube, filter, and instrument can reduce the gas flow.

The device monitors the gas flow in two different modes that are depicted in the status line as follows:

- **Veress insufflation** (1-5 l/min)
- **Insufflation: Standard** (6-40 l/min)

Veress insufflation

Veress insufflation is a gentle type of insufflation that prevents exceeding the preset nominal pressure by the actual pressure even in case of small volumes. To minimize risks in case of a faulty puncture, the manufacturer recommends to begin a procedure (filling abdomen with CO₂) with the setting Veress insufflation.

CAUTION!

Automatic venting system

Please note that the automatic venting system is not active during Veress insufflation mode.

Insufflation: Standard

Start the device by pressing the **Start** function field.

With a nominal flow setting of < 6 l/min, the status line depicts **Veress insufflation**. In case of more than 5 l/min, the display depicts **Insufflation: Standard**.

Stop the device by pressing the **Stop** function field.

NOTE!

Reduced gas flow

Tube, filter, and instrument can reduce the gas flow.

6.5 Performing the Function Test (in Standard Operating Mode)

Check the disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid. Sterilize reusable instruments (**NOT FOR SALE IN USA**) before surgery to avoid infections.

For your own safety and that of your patient, use only original accessories.

WARNING!

Functional test


The functional test must be performed by the user prior to each surgery.

Preparation

Device Check

Fig. 6–1 Layout for the Device Test

- ① ON/OFF button (push-push)
- ② Nominal pressure display
- ③ Start/Stop function field
- ④ Gas consumption display
- ⑤ Nominal flow display
- ⑥ Insufflation tube set
- ⑦ Stopcock (valve)
- ⑧ Veress cannula

1. The device is switched off, no tube set is connected.
2. Make sure the gas supply is connected and open.
3. Use the ON/OFF button (push-push) ① to turn the device on. The device now conducts a device check.
4. Select Standard operating mode in the insufflation mode selection.
5. Connect an original insufflation tube set ⑥ to the device.
6. Attach the insufflation tube to the Veress cannula ⑧.
7. Close the stopcock (valve) ⑦ of the Veress cannula.
8. Press the function field ④ to set the gas consumption display to zero 0 if the display field of the gas consumption depicts a value.
9. Set the nominal pressure to 15 mmHg ② and the gas flow to 3 l/min ⑤.
10. Start insufflation: Press the **Start** ③ function field.
11. The display status line depicts **Occlusion** after max. 4 seconds.
12. Stop insufflation: Press the **Stop** ③ function field.
13. Repeat steps 7 to 11 without Veress cannula and with closed tube end. The previously connected Veress cannula has a leak if gas consumption is now below 0.4 l.

14. Repeat steps 8 to 12 without Veress cannula and without tube if another leak becomes apparent. Close the end of the insufflation tube connection for this test. If the gas consumption is then below 0.4 l, the previously used tube set has a leak.
15. If another leak is detected, this leak is then directly within the device. Make sure the device can no longer be operated until an authorized service technician conducts the appropriate checks.

**WARNING!****Leaky insufflation tube**

Never work with a leaky insufflation tube, accessory, and/or device. This can lead to an incorrect measurement of the actual pressure values, which can cause an uncontrolled pressure increase in the abdomen.

**WARNING!****Leak in the system**

If the actual gas consumption is higher than 0.4 l, there is a leak in the system. If this happens, use steps 13 to 15 outlined below to locate the leak.

Testing gas heating functionality

The gas heating function will be indicated with the gas heating symbol on the display if the "heating tube set" is used.

6.5.1 Filling Tube System with CO₂

At least 1 l of CO₂ must be released from the system before every surgery while the tube is connected and the end of the tube is open to expel any air within the tube system and the device.

1. Start insufflation: Press the **Start** function field. Wait until the gas consumption display shows **1.0 l**.
2. Stop insufflation: Press the **Stop** function field.
3. Press the function field to reset the gas consumption display (see Fig. Screen display [► 30], (13)), to reset the gas consumption display to **0.0 l**. This ensures the gas consumption is correctly displayed during surgery.

The function test is complete. The device is tested and ready to be used for surgery.

**WARNING!****Defect suspected or detected**

Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

6.6 Using the Device during Surgery**WARNING!****Functional test**

The functional test must be successfully completed prior to each surgery.

**CAUTION!****Actual pressure exceeds the set nominal pressure**

The venting system is automatically triggered (not during Veress mode and not if the venting valve has been deactivated in the Configuration menu) if the measured actual pressure exceeds the set nominal pressure during insufflation. The device interrupts insufflation and releases gas until the actual pressure has dropped below the set nominal value.

- The device is turned on.
- An insufflation tube set is connected.

6.6.1 Insufflating with the Veress Cannula

On delivery from the factory, following values are set for Standard operating mode:

- a gas flow value of 3 l/min and
- Veress insufflation for insufflation through the Veress cannula.

CAUTION!

Automatic venting system

Up to 5 l/min Veress insufflation is set automatically. Please note that the automatic venting system is not active during Veress insufflation operating mode (the venting valve is switched off automatically). The gas transport through the device is also specifically customized for use with a Veress cannula.

WARNING!

Position of Veress canula

This instructions for use does not include instructions for the safe use of the Veress cannula. Only after you have verified that a pneumoperitoneum can be generated, should a gas flow of more than 3 l/min and a pressure of more than 10 mmHg be selected. Insert the Veress cannula into the abdomen. Check to see if the Veress cannula is correctly positioned in the abdomen.

1. Attach the insufflation tube (6) (see Fig. Layout for the Device Test [► 39]) to the Veress cannula (8).
2. Select the desired nominal pressure (2) and nominal gas flow (5).
3. Start insufflation: Press the function field **Start** (3). Check the actual pressure display and the gas consumption display.

6.6.2 Insufflating with the Trocar

1. Insert the trocar into the abdomen.
2. Connect the Luer Lock connection of the insufflation tube to the trocar.
3. Make sure the trocar is correctly positioned in the abdomen. Then select the desired pressure and desired gas flow as intraoperative conditions.
4. The actual pressure display shows the current measured value for insufflation. As soon as this value approximates the selected nominal value, the gas flow is automatically minimized. The gas consumption display shows the volume of gas consumed.
5. Check how the patient's body reacts to the selected pressure and gas flow rate. Compare the abdominal filling rate to the selected nominal pressure. You can change the nominal gas flow and the nominal pressure during surgery without interrupting the insufflation process.

6.6.3 Stopping Insufflation

Press the **Stop** function field. The following values are depicted in the display:

- Gas consumption display: last measured value.
- Actual pressure: current measured value
- Actual gas flow: 0 l/min
- Nominal pressure: last set value. In case of exceeding the pressure safety limit, the nominal pressure value will be reset to the pressure safety limit value.
- Nominal gas flow: last set value. In case of exceeding the flow safety limit, the nominal flow value will be reset to the flow safety limit of 5 l/min value (only if activated).
- The display status field depicts **Insufflation stopped**.

en

1. Remove the tube set from the device. Observe applicable hygiene regulations when disposing of the tube set.
 2. Close the gas supply.
 3. Use the ON/OFF button (push-push) to turn the device off.
-

CAUTION!

Leftover fluid

If the tube set remains connected to the device, there is the danger that leftover fluid in the tube or the instruments will be sucked into the device.

NOTE!

Hygiene regulations

Observe applicable hygiene regulations when disposing of the tube set.

7 Device Use and Control with Pediatric Operating Mode

Pediatric operating mode is used to distend the peritoneal cavity during laparoscopic procedures by using a gas supply.

Pediatric operating mode is designed specifically for use on newborns, infants, and children. For use with children, the control of the device - depending on nominal flow selection - is optimized and recommended based on age and weight as per the following:

Age group	Weight	Flow range
Children younger than 1 year	approx. 1-9 kg	0.1-0.5 l/min
Children between 1 and 3 years	approx. 10-15 kg	0.5-1.0 l/min
Children between 3 and 4 years	approx. 16-19 kg	1.0-2.0 l/min
Children between 4 and 14 years	≥ 20 kg	> 2.0 l/min
All children	< 25 kg	≤ 14.0 l/min

Pediatric operating mode is limited to small cavities.

The insufflator limits the pressure to maximal 20 mmHg and the gas flow rate to maximal 20 l/min.

If the nominal flow is set too low, the nominal pressure cannot be reached. Check for possible leaks. Due to the special method used during the **Pediatric** application with the device, the speed of equalizing leaks is slower than when using the High Flow application (lower effective flow during Pediatric application).

Veress insufflation is possible during **Pediatric** operating mode.

The device may not be used to fill an abdomen with CO₂ if a laparoscopy is contraindicated. Please consult the manual of your laparoscope for absolute and relative contraindications. The device is not suitable for hysteroscopic insufflations, i.e., it may not be used to distend the uterus.

The gas flow may not exceed 14 l/min when performing a laparoscopy on infants or patients weighing less than 25 kilos.

Patients

Characteristics

Insufflation parameters

Contraindications

7.1 Device-Specific Risks Associated with Use of the Device in Pediatric Operating Mode

WARNING!

Hypercapnia

Because pediatric patients are especially susceptible to hypercapnia, it is recommended to establish an end-tidal CO₂ monitoring routine.

WARNING!

Procedures with children

Only those who are specially trained and qualified for procedures with children or endoscopic vessel harvesting procedures may use this device for these purposes.

WARNING!

Gas flow limit

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg (approximately 55 US pounds).

WARNING!

Recommended work settings

The flow values listed for laparoscopic procedures performed on newborns, infants, and children are only suggested values. The selection of the suitable flow and pressure values is solely the responsibility of the attending physician.

**WARNING!**

Pneumolabium/pneumoscrotum

Children are at risk of a pneumolabium or pneumoscrotum.

**WARNING!**

Increased airway pressure

When laparoscopic procedures are performed on children, the increased intra-abdominal pressure also increases the risk for higher airway pressures. Always strictly monitor respiration and airway function when performing laparoscopic procedures on children younger than 12 years of age.

**WARNING!**

Compression of the vena cava

When insufflating the abdomen of a child with medical CO₂, an increased risk of compressing the vena cava exists. This risk can be reduced by monitoring the systolic and diastolic blood pressure during the entire surgery.

**WARNING!**


Haemodynamic stability

A laparoscopy performed on children younger than 12 years of age may result in problems of the haemodynamic system due to the CO₂ content in the blood. It is recommended to increase the breathing rate of the patient and to work with low flow values and pressure values not exceeding 12 mmHg. The patient's circulatory system should be monitored at all times.

**WARNING!**

Hypothermia

The insufflation gas flow usually drops significantly after the target pressure has been reached and it is then only required to maintain the abdominal pressure. However, leaks within the abdomen or the instrument can lead to a constant gas flow of above 1 l/min. When operating on children younger than 12, a gas flow of more than 1 l/min poses an increased risk of hypothermia for the patient. Corresponding measures to prevent hypothermia include the use of blankets or pre-warmed gas. The patient's body temperature has to be monitored at all times during surgery.

7.2 Selecting Pediatric Operating Mode

1. Selection of the insufflation mode can be done only while insufflation is stopped. If Pediatric operating mode has not been selected yet, press the function field **Insufflation Mode Selection** for 2 seconds to load the insufflation overview.
2. The display depicts an overview of the insufflation modes. Press the function field **Pediatric** to select the desired insufflation operating mode.
3. The procedure profile is depicted on the display. The displayed parameters correspond with the factory settings or the values set in the configuration menu (see Chapter Configuration Menu (Overview) [► 65]).

7.3 Presetting Nominal Pressure in Pediatric Operating Mode

The nominal pressure can be set during insufflation or while insufflation is stopped. Values may range from 1 to max. 20 mmHg or the value set in the configuration menu of this area.

Increasing/decreasing nominal pressure:

Briefly press the ▲ or ▼ function field to increase or decrease the pressure. Pressing the ▲ or ▼ function field longer than 1.5 seconds activates scrolling in increments of 1.

Select a nominal pressure value between 1 and 20 mmHg. The preset value is indicated in the nominal pressure display.

Pressure safety limit:

Increase the nominal pressure and starting at a value of 12 mmHg the status line depicts **Safety limit** and the nominal pressure value flashes.

The nominal pressure value 12 mmHg is a limit value and should not be exceeded for newborns if at all possible. Pressing the nominal pressure function field ▲ again does not increase the pressure any further.

WARNING!**Exceeding the pressure safety limit**

Exceeding this pressure safety limit is to be decided by and the responsibility of the user/operator.

Release the function field at that point. The display switches back to the nominal value setting after 2 seconds. Now you can set a value up to 15 mmHg. The nominal pressure 15 mmHg is a limit value. This is where the recommended range for the intra-abdominal pressure ends. Pressing the nominal pressure function field ▲ again does not increase the pressure any further. The status line of the display depicts **Safety limit** and the nominal pressure value flashes.

Release the function field at that point. The display switches back to the nominal value setting after 2 seconds. Now you can set a value up to 20 mmHg.

7.4 Presetting Nominal Flow in Pediatric Operating Mode

The nominal flow can be set during insufflation or while insufflation is stopped.

Reducing nominal gas flow:

Press the nominal gas flow function field ▼ to reduce the nominal gas flow. The nominal gas flow is decreased in the range of

- 0.1 l/min to 2 l/min in increments of 0.1 l/min.
- 2 l/min to 20 l/min in increments of 1 l/min.

Pressing the ▼ function field longer than 1.5 seconds activates scrolling.

- In the range of 0.1...2 l/min, rounding down to 1 or 0.1 l/min is applied.
- In the range of 2...20 l/min in increments of 1.

The nominal flow can be set during insufflation or while insufflation is stopped.

Increasing nominal gas flow:

Press the nominal gas flow function field ▲ to increase the nominal gas flow. The nominal gas flow is increased in the range of

- 0.1 l/min to 2 l/min in increments of 0.1 l/min.
- 2 l/min to 20 l/min in increments of 1 l/min.

Pressing the ▲ function field longer than 1.5 seconds activates scrolling.

- In the range of 0.1 l/min...< 2 l/min, rounding up to 1 or 2 l/min is applied.
- Then in one-step increments down to 20 l/min.

Flow safety limit:

When increasing the nominal gas flow, starting with 5 l/min the status line of the display depicts **Safety limit** and the value of the nominal gas flow flashes. (This flow safety limit can be activated or deactivated in the Configuration menu).

Veress insufflation

Release the function field to increase the nominal gas flow further. The display switches back to the nominal value setting after 2 seconds. Now you can set a value of up to 20 l/min.

WARNING!

Exceeding the flow safety limit

Exceeding this flow safety limit is to be decided by and the responsibility of the user/operator.

Veress insufflation operating mode:

The device monitors the gas flow in two different operating modes:

- **Veress insufflation** (≤ 1 l/min)
- **Insufflation: Pediatric** (> 1 l/min)

Veress insufflation is intended for the gentle development of a pneumoperitoneum. The nominal gas flow emitted by the device in this case is very low (in the range of 0.1 l/min to 1 l/min).

To minimize risks in case of a faulty puncture, the manufacturer recommends to begin a procedure (filling abdomen with CO₂) with the setting Veress insufflation.

Start the device by pressing the function field **Start**.

CAUTION!

Automatic venting system

Please note that the automatic venting system is only active during the Veress insufflation mode if "In Veress insufflation ON" has been set in the configuration menu (only possible in Pediatric operating mode, see Chapter Configuration Menu (Overview) [p. 65])

Venting system "In Veress insufflation OFF"

With a nominal flow setting of ≤ 1 l/min, the status line depicts **Veress insufflation** after pressing the Start/Stop function field.

Above 1 l/min, **Insufflation: Pediatric** is displayed.

Venting system "In Veress insufflation ON"

With this setting, pressing the Start function field displays **Insufflation: Pediatric** even if the nominal flow is set to ≤ 1 l/min.

Insufflation: Pediatric

During **Insufflation: Pediatric** any pressure loss due to leaks can be quickly equalized. APC technology (Advanced Pressure Control) enables raising the actual pressure gently to the level of the nominal pressure.

NOTE!

Reduced gas flow

Tube, filter, and instrument can reduce the gas flow.

Preparation

7.5 Performing the Function Test in Pediatric Operating Mode before Using the Device during Surgery

Check the disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid. Sterilize reusable instruments (**NOT FOR SALE IN USA**) before surgery to avoid infections.


For your own safety and that of your patient, use only original accessories.

WARNING!

Functional test

The functional test must be performed by the user prior to each surgery.

Device Check

Fig. 7-1 Layout for the Device Test

- ① ON/OFF button (push-push)
- ② Nominal pressure display
- ③ Start/Stop function field
- ④ Gas consumption display
- ⑤ Nominal flow display
- ⑥ Insufflation tube set
- ⑦ Stopcock (valve)
- ⑧ Veress cannula

1. The device is switched off, no tube set is connected.
2. Make sure the gas supply is connected and open.
3. Use the ON/OFF button (push-push) ① to turn the device on. The device now conducts a device check.
4. Select **Pediatric** operating mode in the insufflation mode selection.
5. Connect an original insufflation tube set ⑥ to the device.
6. Attach the insufflation tube to the Veress cannula ⑧.
7. Close the stopcock (valve) ⑦ of the Veress cannula.
8. Press the function field ④ to set the gas consumption display to zero 0 if the display field of the gas consumption depicts a value.
9. Set the nominal pressure to 8 mmHg ② and the gas flow to 2 l/min ⑤.
10. Start insufflation: Press the **Start** ③ function field.
11. Insufflate for approx. 30 s. The display status line depicts **Occlusion** after maximal 4 seconds.
12. Stop insufflation: Press the **Stop** ③ function field.
13. Repeat steps 7 to 11 without Veress cannula and with closed tube end. The previously connected Veress cannula has a leak if gas consumption is now below 0.4 l.
14. Repeat steps 8 to 12 without Veress cannula and without tube if another leak becomes apparent. Close the end of the insufflation tube connection for this test. If the gas consumption is then below 0.4 l, the previously used tube set has a leak.
15. If another leak is detected, this leak is then directly within the device. Make sure the device can no longer be operated until an authorized service technician conducts the appropriate checks.

**WARNING!****Leaky insufflation tube**

Never work with a leaky insufflation tube, accessory, and/or device. This can lead to an incorrect measurement of the actual pressure values, which can cause an uncontrolled pressure increase in the abdomen.

**WARNING!****Leak in the system**

If the actual gas consumption is higher than 0.4 l, there is a leak in the system. If this happens, use steps 13 to 15 outlined below to locate the leak.

7.5.1 Filling Tube System with CO₂

At least 1 l of CO₂ must be released from the system before every surgery while the tube is connected and the end of the tube is open to expel any air within the tube system and the device.

1. Start insufflation: Press the **Start** function field. Wait until the gas consumption display shows **1.0 l**.
2. Stop insufflation: Press the **Stop** function field.
3. Press the function field to reset the gas consumption display (see Fig. Screen display [► 30], (13)), to reset the gas consumption display to **0.0 l**. This ensures the gas consumption is correctly displayed during surgery.

The function test is complete. The device is tested and ready to be used for surgery.

**WARNING!****Defect suspected or detected**

Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

7.6 Using the Device during Surgery**WARNING!****Functional test**

The functional test must be successfully completed prior to each surgery.

**CAUTION!****Actual pressure exceeds the set nominal pressure**

The venting system is automatically triggered (not during Veress mode and not if the venting valve has been deactivated in the Configuration menu) if the measured actual pressure exceeds the set nominal pressure during insufflation. The device interrupts insufflation and releases gas until the actual pressure has dropped below the set nominal value.

- The device is turned on.
- An insufflation tube set is connected.

7.6.1 Insufflating with the Veress Cannula

On delivery from the factory, following values are set for **Pediatric** operating mode:

- a gas flow value of 0.1 l/min and
- Veress insufflation for insufflation through the Veress cannula.

CAUTION!

Flow value is automatically activated

Please note that with flow values up to 1 l/min the Veress insufflation operating mode is automatically activated (the venting valve is switched off). The gas transport is also specifically customized for use with a Veress cannula.

WARNING!

Safe use of the Veress cannula

These Instructions for use do not include instructions for the safe use of the Veress cannula. Only when you have ensured endoscopically that an aeroperitoneum can be generated, should a gas flow of more than 3 l/min and a pressure of more than 10 mmHg be selected (see corresponding setting for children, Chapter Intended Use [▶ 11]). Insert the Veress cannula into the abdomen. Check to see if the Veress cannula is correctly positioned in the abdomen.

1. Attach the insufflation tube to the Veress cannula.
2. Select the desired nominal pressure and nominal gas flow.
3. Start insufflation: Press the function field **Start**.
Check the actual pressure display and the gas consumption display.

7.6.2 Insufflating with the Trocar

1. Insert the trocar into the abdomen.
2. Connect the Luer Lock connection of the insufflation tube to the trocar.
3. Make sure the trocar is correctly positioned in the abdomen. Then select the desired pressure and desired gas flow as intraoperative conditions.
4. The actual pressure display shows the current measured value for insufflation. As soon as this value approximates the selected nominal value, the gas flow is automatically minimized. The gas consumption display shows the volume of gas consumed.
5. Check how the patient's body reacts to the selected pressure and gas flow rate. Compare the abdominal filling rate to the selected nominal pressure. You can change the nominal gas flow and the nominal pressure during surgery without interrupting the insufflation process.

7.6.3 Stopping Insufflation

Press the **STOP** function field. The following values are depicted in the display:

- Gas consumption display: last measured value.
 - Actual pressure: Current measured value
 - Actual gas flow: 0 l/min
 - Nominal pressure: last set value In case of exceeding the safety limit, the nominal pressure value will be reset to the safety limit value.
 - Nominal gas flow: last set value. In case of exceeding the safety limit, the nominal flow value will be reset to the safety limit value (only if activated).
 - The display status field depicts **Insufflation stopped**.
1. Remove the tube set from the device. Observe applicable hygiene regulations when disposing of the tube set.
 2. Close the gas supply.
 3. Use the ON/OFF switch to turn the device off.

CAUTION!

Leftover fluid

If the tube set remains connected to the device, there is the danger that leftover fluid in the tube or the instruments will be sucked into the device.

en

NOTE!

Hygiene regulations

Observe applicable hygiene regulations when disposing of the tube set.

8 Device Use and Control with Bariatric Operating Mode

Bariatric operating mode is used to distend the peritoneal cavity during laparoscopic procedures by using a gas supply.

Bariatric operating mode is specifically for laparoscopies performed on severely overweight adults (BMI > 30 kg/m²).

Bariatric operating mode enables the rapid insufflation of large volumes.

The insufflator limits the pressure to maximal 30 mmHg and the gas flow rate to maximal 50 l/min.

Veress insufflation is possible during **Bariatric** operating mode.

Patients

Characteristics

Insufflation parameters

8.1 Device-Specific Risks Associated With Use of the Device in Bariatric Operating Mode

WARNING!

Altered Respiratory Physiology

Always monitor the patient's respiratory functions during the entire surgery. The larger body mass supported by the thoracic cage and the larger amount of fat in the abdominal cavity may reduce the elasticity of the thoracic wall. In addition, the increased intra-abdominal pressure secondary to insufflation may alter the normal physiological lung parameters thus resulting in a reduction of the functional lung volume. Shallow, rapid breathing is symptomatic of this condition. Even modest physical stress causes a tremendous increased demand for oxygen, which stands in contrast to the ineffective respiratory musculature that requires more oxygen because it must overcome the reduced elasticity of the thoracic cage. The functional capacity of the lungs is small and even moderate stress can lead to respiratory failure.

WARNING!

Subcutaneous emphysema

When puncturing the thicker abdominal wall of morbidly obese patients with the Veress cannula or the trocar, carefully monitor the correct position of the instrument in the abdomen.

WARNING!

Idiosyncratic reactions

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

WARNING!

CO₂ absorption

CO₂ is absorbed during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

en

WARNING!

CO₂ supersaturation

To avoid generating CO₂ supersaturation, an increased level of respiratory activity is required. An overweight patient's oxygen demand and carbon dioxide production are greater and increase faster under physical stress than do those of patients with normal body weight.

WARNING!

Heart and cardiovascular insufficiency

Constantly monitor all heart and cardiovascular parameters during surgery since morbidly obese patients have an increased risk of heart and cardiovascular insufficiencies.

WARNING!

Metabolic and cardiac reactions

Insufflating CO₂ may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

WARNING!

Hypothermia/monitoring body temperature

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solution not pre-heated to body temperature.

WARNING!

Dehydration

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

WARNING!**Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO₂ embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**WARNING!****Emphysema**

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries (> 200 min.), the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar accesses immediately.


**WARNING!****Additional insufflation sources / automatic venting system**

Make sure the automatic venting system is activated (see Chapter Configuration Menu (Overview) [► 65]). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

8.2 Selecting Bariatric Operating Mode

1. Selection of the insufflation mode can be done only while insufflation is stopped. If Bariatric operating mode is not yet set, press the function field **Insufflation Mode Selection** for 2 seconds (see Chapter Monitor with Touch Screen [► 30]), to load the insufflation overview.
2. The display depicts an overview of the insufflation modes. Select the desired insufflation mode by pressing the **Bariatric** function field.
3. The procedure profile is depicted on the display. The displayed parameters correspond with the factory settings or the values set in the configuration menu (see Chapter Configuration Menu I).

8.3 Presetting Nominal Pressure in Bariatric Operating Mode

The nominal pressure can be set during insufflation or while insufflation is stopped. Values may range of 1 to max. 30 mmHg or correspond with the value set in the configuration menu in this area.

Increasing/decreasing nominal pressure:

Briefly press the ▲ or ▼ function field to increase or decrease the pressure. Pressing the ▲ or ▼ function field depressed longer than 1.5 seconds activates scrolling in increments of 1.

Pressure safety limit:

When **increasing** the nominal pressure, the status line of the display depicts **Safety Limit** starting at 15 mmHg and the nominal pressure value flashes. The nominal pressure value of 15 mmHg is a pressure safety limit value. This is where the recommended range for the intra-abdominal pressure ends. Pressing the nominal pressure function field ▲ again does not increase the pressure any further. Release the function field at that point. The display switches back to the nominal value setting after 2 seconds. Now you can set a value up to 30 mmHg.

**CAUTION!****Exceeding pressure safety limit**

Exceeding this pressure safety limit is to be decided by and the responsibility of the user/ operator.

8.4 Presetting Nominal Flow in Bariatric Operating Mode

The nominal flow can be increased or decreased during insufflation or while insufflation is stopped.

1. Briefly press the ▲ or ▼ function field to increase or decrease the flow. The nominal flow can be increased from 1 to max. 50 l/min.
 - Briefly press the corresponding field to set values in increments of 1.
 - Pressing the ▲ or ▼ function field depressed longer than 1.5 seconds activates scrolling through the possible gas flow values.
2. The preset values of the configuration menu can be changed individually (see Chapter Configuration Menu I [► 66]). Select a nominal gas flow between 1 and 50 l/min. The preset value is displayed in the nominal flow display. The values for the nominal gas flow refer to a device without connected tube, filter, or instrument. Tube, filter, and instrument can reduce the gas flow.

The device monitors the gas flow in two different modes:

- **Veress insufflation (1 to 5 l/min)**
- **Insufflation: Bariatric (6 - 50 l/min)**

Veress insufflation is a gentle type of insufflation that prevents exceeding the preset nominal pressure by the actual pressure even in case of small volumes. To minimize risks in case of a faulty puncture, the manufacturer recommends to begin a procedure (filling abdomen with CO₂) with the setting Veress insufflation.

**CAUTION!****Automatic venting system**

Please note that the automatic venting system is not active during Veress insufflation mode.

Insufflation: Bariatric

With **Insufflation: Bariatric**, any pressure loss due to leaks can be quickly equalized. APC technology (Advanced Pressure Control) enables raising the actual pressure gently to the level of the nominal pressure. In case of large volumes, the actual pressure usually does not exceed the nominal pressure (see Chapter Safety Functions [► 72]).

Start the device by pressing the **Start** function field.

Veress Insufflation is displayed with a nominal flow setting of < 6 l/min. Above 5 l/min **Insufflation: Bariatric** is displayed.

Stop the device by pressing the **Stop** function field.

**NOTE!****Reduced gas flow**

Tube, filter, and instrument can reduce the gas flow.

8.5 Performing the Function Test in Bariatric Operating Mode before Using the Device during Surgery

Check the disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid. Sterilize reusable instruments (**NOT FOR SALE IN USA**) before surgery to avoid infections.

For your own safety and that of your patient, use only original accessories.

Preparations

WARNING!


Functional test

The functional test must be performed by the user prior to each surgery.

Device Check

Fig. 8–1 Layout for the Device Test

- ① ON/OFF button (push-push)
- ② Nominal pressure display
- ③ Start/Stop function field
- ④ Gas consumption display
- ⑤ Nominal flow display
- ⑥ Insufflation tube set
- ⑦ Stopcock (valve)
- ⑧ Veress cannula

1. The device is switched off, no tube set is connected.
2. Make sure the gas supply is connected and open.
3. Use the ON/OFF button (push-push) ① to turn the device on. The device now conducts a device check.
4. Select Bariatric mode in the insufflation mode selection.
5. Connect an original insufflation tube set ⑥ to the device.
6. Attach the insufflation tube to the Veress cannula ⑧.
7. Close the stopcock (valve) ⑦ of the Veress cannula.
8. Press the function field ④ to set the gas consumption display to zero 0 if the display field of the gas consumption depicts a value.
9. Set the nominal pressure to 15 mmHg ② and the gas flow to 5 l/min ⑤.
10. Start insufflation: Press the **Start** ③ function field.
11. Insufflate for approx. 30 s. The display status line depicts Occlusion after max. 4 seconds.
12. Stop insufflation: Press the **Stop** ③ function field.
13. Repeat steps 7 to 11 without Veress cannula and with closed tube end. The previously connected Veress cannula has a leak if gas consumption is now below 0.4 l.
14. Repeat steps 8 to 12 without Veress cannula and without tube if another leak becomes apparent. Close the end of the insufflation tube connection for this test. If the gas consumption is then below 0.4 l, the previously used tube set has a leak.
15. If another leak is detected, this leak is then directly within the device. Make sure the device can no longer be operated until an authorized service technician conducts the appropriate checks.

WARNING!

Leaky insufflation tube

Never work with a leaky insufflation tube, accessory, and/or device. This can lead to an incorrect measurement of the actual pressure values, which can cause an uncontrolled pressure increase in the abdomen.

WARNING!

Leak in the system

If the actual gas consumption is higher than 0.4 l, there is a leak in the system. If this happens, use steps 13 to 15 outlined below to locate the leak.

Testing gas heating functionality

The gas heating function will be indicated with the gas heating symbol on the display if the "heating tube set" is used.

8.5.1 Filling Tube System with CO₂

At least 1 l of CO₂ must be released from the system before every surgery while the tube is connected and the end of the tube is open to expel any air within the tube system and the device.

1. Start insufflation: Press the **Start** function field. Wait until the gas consumption display shows **1.0 l**.
2. Stop insufflation: Press the **Stop** function field.
3. Press the function field to reset the gas consumption display (see Fig. Screen display [► 30], (13)), to reset the gas consumption display to **0.0 l**. This ensures the gas consumption is correctly displayed during surgery.

The function test is complete. The device is tested and ready to be used for surgery.

WARNING!

Defect suspected or detected

Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

8.6 Using the Device during Surgery

WARNING!

Functional test

The functional test must be successfully completed prior to each surgery.

CAUTION!

Actual pressure exceeds the set nominal pressure

The venting system is automatically triggered (not during Veress mode and not if the venting valve has been deactivated in the Configuration menu) if the measured actual pressure exceeds the set nominal pressure during insufflation. The device interrupts insufflation and releases gas until the actual pressure has dropped below the set nominal value.

- The device is turned on.
- An insufflation tube set is connected.

8.6.1 Insufflating with the Veress Cannula

On delivery from the factory, following values are set for **Bariatric** operating mode:

- a gas flow value of 5 l/min and
- Veress insufflation for insufflation through the Veress cannula.

CAUTION!**Automatic venting system**

Up to 5 l/min Veress insufflation is set automatically. Please note that the automatic venting system is not active during Veress insufflation operating mode (the venting valve is switched off automatically). The gas transport through the device is also specifically customized for use with a Veress cannula.

**WARNING!****Position of Veress cannula**

This instructions for use does not include instructions for the safe use of the Veress cannula. Only after you have verified that a pneumoperitoneum can be generated, should a gas flow of more than 3 l/min and a pressure of more than 10 mmHg be selected. Insert the Veress cannula into the abdomen. Check to see if the Veress cannula is correctly positioned in the abdomen.

1. Attach the insufflation tube (6) (see Fig. Layout for the Device Test) the Veress cannula (8).
2. Select the desired nominal pressure (2) and nominal gas flow (5).
3. Start insufflation: Press the function field **Start** (3).

Check the actual pressure display and the gas consumption display.

8.6.2 Insufflating with the Trocar

1. Insert the trocar into the abdomen.
2. Connect the Luer Lock connection of the insufflation tube to the trocar.
3. Make sure the trocar is correctly positioned in the abdomen. Then select the desired pressure and desired gas flow as intraoperative conditions.
4. The actual pressure display shows the current measured value for insufflation. As soon as this value approximates the selected nominal value, the gas flow is automatically minimized. The gas consumption display shows the volume of gas consumed.
5. Check how the patient's body reacts to the selected pressure and gas flow rate. Compare the abdominal filling rate to the selected nominal pressure. You can change the nominal gas flow and the nominal pressure during surgery without interrupting the insufflation process. The Luer Lock connection of the measuring tube can be connected to another trocar that also inserted when using a tube set with measuring line. This enables the continuous measuring of the pressure.

8.6.3 Stopping Insufflation

Press the **Stop** function field. The following values are depicted in the display:

- Gas consumption display: last measured value.
- Actual pressure: Current measured value
- Actual gas flow: 0 l/min
- Nominal pressure: last set value In case of exceeding the pressure safety limit, the nominal pressure value will be reset to the pressure safety limit value.
- Nominal gas flow: last set value.

The display status field depicts **Insufflation stopped**.

1. Remove the tube set from the device. Observe applicable hygiene regulations when disposing of the tube set.
2. Close the gas supply.
3. Use the ON/OFF switch to turn the device off.

en

CAUTION!

Leftover fluid

If the tube set remains connected to the device, there is the danger that leftover fluid in the tube or the instruments will be sucked into the device.

NOTE!

Hygiene regulations

Observe applicable hygiene regulations when disposing of the tube set.

9 Device Use and Control in Vessel Harvest Operating Mode

Vessel Harvest mode is used to distend the cavity along the vena saphena magna and radial artery during an endoscopic vessel harvesting procedure. Carrying out the surgery for vessel harvesting requires the use of a special instrument.

Vessel Harvest mode is suitable for patients of all ages regardless of weight.

Vessel Harvest mode is limited to small cavities.

The insufflator limits the pressure to maximal 20 mmHg and the gas flow rate to maximal 10 l/min.

Veress insufflation is not possible during **Vessel Harvest** mode.

The device may not be used for the endoscopic harvesting of vessels if this surgical application is contraindicated. Please consult the manual of that instrument for additional information and special application indications.

Patients

Characteristics

Insufflation parameters

Contraindications

9.1 Device-Specific Risks Associated with Use of the Device in Vessel Harvest Operating Mode

WARNING!

Procedures with children

Only those who are specially trained and qualified for procedures with children or endoscopic vessel harvesting procedures may use this device for these purposes.

WARNING!

Instrument used for CO₂ insufflation

Before using the insufflator to endoscopically harvest vessels, please check whether the instrument used is intended for CO₂ insufflation.

WARNING!

Pneumoperitoneum

When a vessel is harvested from the leg of a patient with a perforated groin, it is possible for CO₂ to reach the abdomen and cause a pneumoperitoneum. Make sure the abdomen does not fill with CO₂ during surgery.

WARNING!

Idiosyncratic reactions

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

WARNING!

CO₂ absorption

Due to the special surgical procedures - start of the heart bypass operation, and the endoscopic removal of the vessel - special care has to be taken as CO₂ is always absorbed through the tissue of the patient during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption.

**WARNING!****Metabolic and cardiac reactions**

Due to the special surgical conditions - start of the heart bypass surgery and vessel harvesting - it is especially important to remember the existing risk of metabolic acidosis when insufflating with CO₂. This can lead to cardiac irregularities expressed with the following symptoms:


- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO₂ embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**9.2 Selecting Vessel Harvest Operating Mode**

1. Selection of the insufflation mode can be done only while insufflation is stopped. If Endoscopic Vessel Harvesting Operating Mode is not yet set, press the function field **Insufflation Mode Selection** for 2 seconds to load the insufflation overview.
2. The display depicts an overview of the insufflation modes. Press the function field **Vessel Harvest** (endoscopic procedure) to select the desired insufflation operating mode.
3. The procedure profile is depicted on the display. The displayed parameters correspond with the factory settings or the values set in the configuration menu (see Chapter Configuration Menu I [► 66]).

9.3 Presetting Nominal Pressure in Vessel Harvest Operating Mode

The nominal pressure can be set during insufflation or while insufflation is stopped. Values may be in the range of 1 to max. 20 mmHg or correspond to the value set in the configuration menu in this area.

Increasing/decreasing nominal pressure:

Briefly press the ▲ or ▼ function field to increase or decrease the pressure. Pressing the ▲ or ▼ function field longer than 1.5 seconds activates scrolling in increments of 1.

Pressure safety limit:

If you increase the nominal pressure and starting at a value of 12 mmHg, the status line depicts **Safety limit** and the nominal pressure value flashes.

The nominal pressure value 12 mmHg is a limit value and should not be exceeded when performing endoscopic vessel harvesting surgery. Pressing the nominal pressure function field ▲ again does not increase the pressure any further.

WARNING!

Exceeding the pressure safety limit

Exceeding this pressure safety limit is to be decided by and the responsibility of the user/operator.

Release the function field at that point. The display switches back to the nominal value setting after 2 seconds. Now you can set a value up to 15 mmHg.

Starting at 15 mmHg, the status line again depicts **Safety limit** and the nominal pressure value flashes.

The nominal pressure value 15 mmHg is another limit value and should only be exceeded when performing vessel harvesting surgery if special circumstances require doing so. Pressing the nominal pressure function field ▲ again does not increase the pressure any further.

Release the function field at that point. The display switches back to the nominal value setting after 2 seconds. Now you can set a value up to 20 mmHg.

9.4 Presetting Nominal Flow in Vessel Harvest Operating Mode

The nominal flow can be set during insufflation or while insufflation is stopped.

Increasing nominal gas flow:

- Press the function field ▲ to increase the flow. You can increase the nominal flow from **1 to maximal 10 l/min**.
- Briefly press the corresponding function field to set values in increments of 0.5 for a range from 1 to 5 l/min.
- Set values in increments of 1 for a range from 5 to 10 l/min.
- Pressing the ▲ or ▼ function field longer than 1.5 seconds activates scrolling through the gas flow rates 1, 4, 10 l/min (in the Configuration menu).

Reducing nominal gas flow:

Briefly press the nominal gas flow function field ▼ to reduce the nominal gas flow.

Flow safety limit:

The flow safety limit can be activated or deactivated in the Configuration menu.

Increasing the nominal flow, and starting at a value of 6 l/min the status line depicts **Safety limit** and the nominal flow value flashes.

Release the function field to increase the nominal gas flow further. The display switches back to the nominal value setting after 2 seconds. Now you can set a value of up to 10 l/min.

WARNING!

Exceeding the flow safety limit

Exceeding this flow safety limit is to be decided by and the responsibility of the user/operator.

9.5 Performing the Function Test in Vessel Harvest Operating Mode before Using the Device during Surgery

Check the disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid. Sterilize reusable instruments (**NOT FOR SALE IN USA**) before surgery to avoid infections.

For your own safety and that of your patient, use only original accessories.


Preparation

Device check

Fig. 9–1 Layout for the Device Test

- ① ON/OFF button (push-push)
- ② Nominal pressure display
- ③ Start/Stop function field
- ④ Nominal flow display
- ⑤ Insufflation tube set
- ⑥ Luer Lock connection
- ⑦ Gas consumption display

1. The device is switched off, no tube set is connected.
2. Make sure the gas supply is connected and open.
3. Use the ON/OFF button (push-push) ① to turn the device on. The device now conducts a device check.
4. Select **Vessel Harvest** operating mode in the insufflation mode selection.
5. Connect an original insufflation tube set ⑤ to the device.
6. Press the function field for the gas consumption display ⑦ to set the display to **0 l**.
7. Select a nominal pressure of 10 mmHg ② and a gas flow of 1 l/min ④.
8. Close off the end of the tube (Luer Lock) ⑥.
9. Start insufflation: Press the **Start** ③ function field.
10. Insufflate for approx. 30 s. The display status line depicts Occlusion after maximal 4 seconds.
11. Stop insufflation: Press the **Stop** ③ function field.
12. Replace the tube set and repeat the test from step 6 to step 11.
13. If another leak is detected, this leak is then directly within the device. Make sure the device can no longer be operated until an authorized service technician conducts the appropriate checks.

WARNING!

Leak in the system

If the actual gas consumption is higher than 0.4 l, there is a leak in the system.

WARNING!

Leaky insufflation tube

Never work with a leaky insufflation tube, accessory, and/or device. This can lead to an incorrect measurement of the actual pressure values, which can cause an uncontrolled pressure increase in the abdomen.

9.5.1 Filling Tube System with CO₂

At least 1 l of CO₂ must be released from the system before every surgery while the tube is connected and the end of the tube is open to expel any air within the tube system and the device.

1. Start insufflation: Press the **Start** function field. Wait until the gas consumption display shows **1.0 l**.
2. Stop insufflation: Press the **Stop** function field.
3. Press the function field to reset the gas consumption display (see Fig. Screen display [► 30], (13)), to reset the gas consumption display to **0.0 l**. This ensures the gas consumption is correctly displayed during surgery.

The function test is complete. The device is tested and ready to be used for surgery.

WARNING!

Defect suspected or detected

Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

9.6 Using the Device during Surgery

WARNING!

Functional test

The functional test must be performed by the user prior to each surgery.

CAUTION!

Actual pressure exceeds the set nominal pressure

The venting system is automatically triggered (only if activated in the Configuration menu) if the measured actual pressure exceeds the set nominal pressure during insufflation. The device interrupts insufflation and releases gas until the actual pressure has dropped below the set nominal value.

- The device is turned on.
- An insufflation tube set is connected.

9.6.1 Insufflation with Endoscopic Vessel Harvest Instrument

Set the **Vessel Harvest** operating mode in the configuration menu for an endoscopic vessel harvesting procedure. If the device was switched off during a vessel harvesting procedure, the display depicts **Vessel Harvest** after restarting the device indicating that the device is already using the Vessel Harvest mode.

CAUTION!

Instrument in conjunction with CO₂

Before you use the insufflator for a vessel harvesting procedure, make sure to read the manual of your endoscopic instrument to check whether the instrument can be used in conjunction with CO₂.

1. Connect the special vessel harvesting instrument with the insufflation tube (5) (see Fig. Layout for the Device Test [► 62]).
Use a manufacturer's recommended instrument to harvest vessels.
2. Select the necessary nominal pressure (2) and nominal gas flow (4).
3. Connect the Luer lock (6) connection of the insufflation tube to the vessel harvesting instrument.
4. Select the desired pressure and desired gas flow as intraoperative conditions.

5. Start insufflation: Press the **Start** function field (3). Check the actual pressure display and the gas consumption display.
6. The actual pressure display shows the current measured value for insufflation. As soon as this value approximates the selected nominal value, the gas flow is automatically minimized. The gas consumption display shows the volume of gas consumed.
7. Check how the patient's body reacts to the selected pressure and gas flow rate. You can change the nominal gas flow and the nominal pressure during surgery without interrupting the insufflation process.

9.6.2 Stopping Insufflation

Press the function field **Stop**. The display depicts the following values:

- Gas consumption display: last measured value.
- Actual pressure: Current measured value
- Actual gas flow: 0 l/min
- Nominal pressure: last set value. In case of exceeding the pressure safety limit, the nominal pressure value will be reset to the pressure safety limit value.
- Nominal gas flow: last set value. In case of exceeding the flow safety limit, the nominal flow value will be reset to the flow safety limit value (only if activated).

The status field of the display depicts **Insufflation stopped** followed by **Mode: Vessel Harvest** alternately.

1. Remove the tube set from the device. Observe applicable hygiene regulations when disposing of the tube set.
2. Close the gas supply.
3. Use the ON/OFF switch to turn the device off.

CAUTION!

Leftover fluid

If the tube set remains connected to the device, there is the danger that leftover fluid in the tube or the instruments will be sucked into the device.

NOTE!

Hygiene regulations

Observe applicable hygiene regulations when disposing of the tube set.

10 Configuration Menu (Overview)

Device parameters are changed with the configuration menu. The following is an overview and details on the following pages (factory settings = **bold**).

Configuration menu

Configuration Menu I				
Menu level	↔	1. Submenu level	2. Submenu level	3. Submenu level
First nominal pressure ↑↓		Standard: 15 mmHg (Setting range: 1 to 15 mmHg) Bariatric: 15 mmHg (Setting range: 1 to 15 mmHg) Vessel Harvest: 10 mmHg (Setting range 1-12 mmHg) Pediatric: 8 mmHg (Setting range 1-12 mmHg)		
Venting system ↑↓	↔	Venting valve status	↔	Standard:* With Veress insufflation OFF Venting system OFF
			↔	Bariatric:* With Veress insufflation OFF Venting system OFF
			↔	Vessel Harvest:* Venting system ON Venting system OFF
			↔	Pediatric:* With Veress insufflation ON With Veress insufflation OFF Venting system OFF
	↔	Venting valve pressure limit	↔	2 mmHg, 3 mmHg , 4 mmHg, 5 mmHg
	↔	Venting valve response time	↔	2 s, 3 s , 4 s, 5 s
Gas supply ↑↓	↔	House gas Gas bottle		
Acoustic signal volume ↑↓	↔	Level 1 Level 2 Level 3		

*On the insufflator screen, inactive menu items are displayed in gray.

Configuration menu II				
Menu level	↔	1. Submenu level	2. Submenu level	3. Submenu level
Gas flow rate* ↑↓ ↑↓	↔	Standard: Level 1 = 3 l/min (Setting range: 1-20 l/min) Level 2 = 20 l/min (Setting range: Level 1, max. 40 l/min) Level 3 = 40 l/min (Setting range: Level 2, max. 40 l/min)		
		Bariatric: Level 1 = 5 l/min (Setting range: 1-20 l/min) Level 2 = 25 l/min (Setting range: Level 1, max. 50 l/min) Level 3 = 50 l/min (Setting range: Level 2, max. 50 l/min)		
	↔	Vessel Harvest: Level 1 = 1 l/min (range 1-10) Level 2 = 4 l/min (range level 1-10) Level 3 = 10 l/min (range level 2-10)		

Configuration menu II			
		Pediatric: inactive	
First nominal gas flow	↔	Pediatric: 0.1 l/min (Setting range 0.1-1 l/min) - active only in Pediatric operating mode	
Maximum nominal pressure ↑↓	↔	Standard: 30 mmHg (Setting range 5-30 mmHg) Bariatric: 30 mmHg (Setting range 5-30 mmHg) Vessel Harvest: 20 mmHg (Setting range 5-20 mmHg) Pediatric: 20 mmHg (Setting range 5-20 mmHg)	
Flow safety limit ↑↓	↔	Limit ON Limit OFF Active only in Pediatric operating mode and Endoscopic Vessel Harvest operating mode	
Occlusion warning signal ↑↓	↔	Standard, Bariatric, Pediatric: Signal ON , Signal OFF Vessel Harvest: Signal ON, Signal OFF	
*On the insufflator screen, inactive menu items are displayed in gray.			

Utility menu						
Menu level	↔	1. Submenu level		2. Submenu level		3. Submenu level
Brightness		100 %, 75 %, 50 %, 25 %				
↑↓						
Language	↔	English Français Deutsch Español Português Italiano	↔	Nederlands Norsk Suomi Greek Svenska Dansk	↔	Polski Română Simplified Chinese Korean Japanese
↑↓						
Software version	↔	Shows current software version				
↑↓						
Service		Access only for service technicians				
*On the insufflator screen, inactive menu items are displayed in gray.						

10.1 Configuration Menu I

Functions of the configuration menu

While insufflation is stopped, press the **Menu** function field (see Fig. Screen display [► 30], (4)) for 2 seconds to open Configuration Menu I.

Depiction/functions in Configuration menu	
	Press the arrow forward function field to access the next menu on the same level.
	Press the arrow backward function field to access the previous menu on the same level without saving settings or to access the previous menu of the next higher menu level.
	Press the (SAVE) function field to save settings. The display depicts (SAVED) for 2 seconds. After saving, you return automatically to the previous level.
	Press the (EXIT) function field to exit the menu or return to the work screen.

10.1.1 Setting First Nominal Pressure

Mode	Default setting	Setting range
Bariatric	15 mmHg	1-15 mmHg
Standard	15 mmHg	1-15 mmHg
Vessel Harvest	10 mmHg	1-12 mmHg
Pediatric	8 mmHg	1-12 mmHg

In Configuration Menu I, press the **First Nominal Pressure** function field to access the setting.

The display switches to the **First Nominal Pressure** menu.

Press the ▲ or ▼ function field to set the **First Nominal pressure**.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the ⇐ function field to return to the previous menu level without saving.
3. Or press **(EXIT)** to return to the work screen without saving.

10.1.2 Setting the Venting Valve

In Configuration Menu I, press the **Venting Valve** function field to access the venting system selection.

The display switches to the **Venting controls** selection menu.

Default settings

Mode	Default setting	Setting range
Bariatric	With Veress insufflation OFF	With Veress insufflation OFF, Venting system OFF
Standard	With Veress insufflation OFF	With Veress insufflation OFF, Venting system OFF
Vessel Harvest	Venting system OFF	Venting system ON, Venting system OFF
Pediatric	With Veress insufflation ON	With Veress insufflation ON, with Veress insufflation OFF, Venting system OFF

Setting the venting valve status

For example, press the **Venting Valve Status** function field.

The display switches to the **Venting Valve Status** selection menu.

For example, press the **Venting System OFF** function field.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the ⇐ function field to return to the previous menu level without saving.
3. Or press **(EXIT)** to return to the work screen without saving.

Default setting = 3 mmHg (for all operating modes)

For example, press the **Venting Valve Pressure Limit** function field.

The display switches to the **Pressure Limit** selection menu.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.

Setting the venting valve pressure limit

Setting the venting system response time

2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

Default setting = 3 s (for all operating modes)

For example, press the **Venting Valve Response Time** function field.

The display switches to the **Venting Valve Response Time** selection menu.

For example, press the **4 s** function field.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.1.3 Setting the Gas Supply Type

Use this menu to select the type of connected gas supply.

Select **House Gas** if you are working with a house gas supply (use the corresponding house gas device adapter available as an additional accessory). The gas supply indicators and displays are described in Chapter Safety Functions [► 72]. If house gas supply has been selected but bottled gas (pressure >15 bar) is connected to the device, the gas supply display automatically switches accordingly.

Select **Bottle** if you want to work with a gas cylinder/bottle. The gas supply indicators and displays are described in Chapter Safety Functions [► 72]. It is not possible to operate the device if **Bottle** gas is selected and a house gas supply is actually connected.

Default setting = House gas

In Configuration Menu I, press the **Gas Supply** function field to access the gas supply selection menu.

The display switches to the **Gas supply** selection menu.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.1.4 Setting the Acoustic Signal Volume

Default setting = Level 3

In Configuration Menu I, press the **Acoustic signal volume** function field to access the volume selection.

The display switches to the **Acoustic signal volume** selection menu. The selection/setting applies to all operating modes.

For example, press the **Level 2** function field to lower the volume.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.2 Configuration Menu II

While insufflation is stopped, press the **Actual** function field for 2 seconds to access the Configuration menu. In the Configuration menu I, press the $\Rightarrow$ to access the Configuration menu II.

Configuration menu II

10.2.1 Setting the Gas Flow Rate

Level	Mode	Default setting	Setting range
Level 1	Bariatric	5 l/min	1-20 l/min
	Standard	3 l/min	1-20 l/min
	Vessel Harvest	1 l/min	1-10 l/min
	Pediatric	Not available	
Level 2	Bariatric	25 l/min	Level 1-50 l/min
	Standard	20 l/min	Level 1-40 l/min
	Vessel Harvest	4 l/min	Level 1-10 l/min
	Pediatric	Not available	
Level 3	Bariatric	50 l/min	Level 2-50 l/min
	Standard	40 l/min	Level 2-40 l/min
	Vessel Harvest	10 l/min	Level 2-10 l/min
	Pediatric	Not available	

In Configuration Menu II, press the **Gas Flow Rates** function field to access the selection menu.

The display switches to the **Gas Supply** selection menu.

For example, press the function field (**flow rate 1: 3 l/min**) to set gas flow rate 1.

The display switches to the Setting menu for gas flow rate 1.

Press the function field $\blacktriangle$ or $\blacktriangledown$ to set the value.

You can now do the following:

1. Save the selected value by pressing the (**SAVE**) function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving.
3. Or press (**EXIT**) to return to the work screen without saving.

10.2.2 Setting the Maximum Nominal Pressure

Mode	Default setting	Setting range
Bariatric	30 mmHg	5-30 mmHg
Standard	30 mmHg	5-30 mmHg
Vessel Harvest	20 mmHg	5-20 mmHg
Pediatric	20 mmHg	5-20 mmHg

In Configuration Menu I, press $\Rightarrow$ to access Configuration Menu II.

The display changes to the configuration menu II.

In Configuration Menu II, press the **Maximal Nominal Pressure** function field.

The display switches to the **Maximal Nominal Pressure** selection menu.

Press the $\blacktriangle$ or $\blacktriangledown$ function field to set the value.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.2.3 Setting the Occlusion Warning Signal

Mode	Default setting
Bariatric	Signal ON
Standard	Signal ON
Vessel Harvest	Signal OFF
Pediatric	Signal ON

In Configuration Menu I, press $\Rightarrow$ to access Configuration Menu II.

The display changes to the configuration menu II.

In the configuration menu II, press the **Occlusion Warning Signal** function field to access the selection menu.

The display switches to the selection menu for the **Occlusion Warning Signal**.

For example, press the **Signal OFF** function field.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.3 Settings Menu

While insufflation is stopped, press the **Actual** function field for 2 seconds to access the Configuration menu. In the Configuration menu I, press the $\Rightarrow$ function field to access the Configuration menu II and there press the $\Rightarrow$ function field again to access the **Utility menu**.

10.3.1 Checking Software Version

In the Settings menu, press the **Software Version** function field to open the corresponding screen.

The display depicts the **Software Version**.

You can now do the following:

1. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
2. Or press **(EXIT)** to return to the work screen without saving.

10.3.2 Brightness

Default setting = 75%

In the Settings menu, press the **Brightness** function field to access the selection menu.

The display switches to the **Brightness** selection menu.

Press the appropriate function field to adjust the brightness (25%, 50%, 75%, or 100%).

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.3.3 Select Language

While insufflation is stopped, press the **Actual** function field for 2 seconds to access the Configuration menu. In Configuration Menu I, press the $\Rightarrow$ function field to access Configuration Menu II and there press the $\Rightarrow$ function field again to access the Setting menu.

Default setting = English

In the Setting menu, press the **Language** function field to access the language selection screen.

The following languages are available:

1. Submenu level	2. Submenu level	3. Submenu level
English	Nederlands	Polski
Français	Norsk	Română
Deutsch	Suomi	Simplified Chinese
Español	Greek	Korean
Português	Svenska	Japanese
Italiano	Dansk	

Press the **English** function field to select English as the display language, for example.

You can now do the following:

1. Save the selected value by pressing the **(SAVE)** function field. The display depicts **(SAVED)** for 2 seconds. After saving, you return automatically to the previous menu level.
2. Press the $\Leftarrow$ function field to return to the previous menu level without saving,
3. Or press **(EXIT)** to return to the work screen without saving.

10.3.4 Service Menu

Access to the service menu is restricted for trained and authorized service personnel. Please consult the service manual for further information.

Automatic venting system

Insufflation source

Exceeding nominal pressure

Exceeding nominal pressure for more than 5 seconds

Nominal pressure limit 30 mmHg

With deactivated venting system or during Veress insufflation

Warning message "Occlusion"

Error message "Contamination"

Contamination/Call service

11 Safety Functions

The device is equipped with an automatic venting system.

When the insufflator detects an overpressure it automatically activates the venting system. The venting system releases gas from the abdomen or from the vessel until the set nominal value has been reached again.

The function "Automatic venting system" can be activated/deactivated in the Configuration menu.

CAUTION!**Automatic venting system**

Please note that the automatic venting system is only active during the Veress insufflation mode if "With Veress Insufflation ON" has been set in the configuration menu.

WARNING!**Limited venting system**

The venting rate of the automatic venting system is limited. Always monitor the actual pressure when using additional insufflation sources.

The manufacturer advises against using additional, non-pressure controlled insufflation sources during minimally invasive surgical procedures.

The use of CO₂-cooled lasers and argon beamers can lead to values exceeding the recommended and safe pressure rating.

After 2 to 5 seconds (depending on the venting valve response time setting in the Configuration menu), the display depicts **Venting System Active** when the nominal pressure has been exceeded by 2-5 mmHg (depending on the venting valve response time setting in the Configuration menu).

If the overpressure cannot be reduced by the automatic venting system within 5 seconds, the display depicts **Overpressure** followed by **Venting System Active**. An acoustic warning is emitted.

When reaching/exceeding the nominal pressure threshold (Standard/Bariatric -> 30 mmHg), the display depicts **Overpressure**. An acoustic warning signal is emitted at the same time.

If the actual pressure exceeds the nominal pressure by 4 mmHg longer than 3 seconds, the display depicts **Overpressure**.

When the tube, Veress cannula, or trocar experiences a temporary blockage, Occlusion is depicted. An acoustic warning signal is emitted. The actual pressure and actual flow displays depict **0**.

The acoustic warning signal can be activated/deactivated in the Configuration menu.

When fluid has penetrated the device via the insufflation tube connection, the display depicts **Contamination** and 3 acoustic warning signals are emitted.

Message is repeated with each Start/Stop.

It is possible to conclude any currently ongoing surgery with this device. Insufflation is not possible after turning the device off and back on using the ON/OFF button (push-push). This is to prevent cross-contamination.

The display depicts **Contamination** alternating with Call service if you are switching on an already contaminated device. The device can no longer be used. The contaminated device has to be clearly marked as contaminated and sealed in two separate protective layers of safety foil. Make sure the device can no longer be used until a qualified service technician conducts the appropriate tests and repairs.

CAUTION!**Venting valve**

The venting valve is no longer available to reduce high pressures after a problem with the venting valve has been detected.

Three beeps and the message **Venting valve defective** are repeated with each Start/Stop.

It is still possible to use the device but pay special attention to the pressure values while this error persists. A service technician should check/repair the device after the completed surgery.

If the power fails for less than 1 second, all settings are retained. If the power fails for an extended time period, the device will function as it does when it is being re-started.

If a device malfunction occurs that prohibits any further use of the device or diminished its safety, the display depicts only

Text message - >Call service

The following text messages may be displayed:

- Contamination
- Electronic defective
- Sensor defective
- Venting valve defective
- Calibration error
- Temperature error

Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

Error message when starting with defective venting valve

Power failure

12 Care and Maintenance

Special care is necessary when servicing, maintaining, and storing the device and its accessories to maintain the functionality of the device and its accessories.

NOTE!

Service or maintenance work

Service or maintenance work may not be carried out during surgery.

NOTE!

Phenol- or chlorine-based disinfectants

Do not use phenol- or chlorine-based disinfectants or cleaning solutions.

CAUTION!

Sterilization of the high pressure hoses

High pressure hoses cannot be sterilized!

Clean the surface of the gas supply hoses at least once a month by wiping with a surface disinfectant based on the active ingredient of aldehydes, quaternary ammonium compounds, or alcohol.

Manufacturer's specifications

12.2 Cleaning the Device

The manufacturer recommends that the device is cleaned after every procedure as follows:

1. Use the **ON/OFF** button (push-push) to turn the device off.
2. Remove the power cable.
3. Wipe the surface of the device with a soft lint-free cloth moistened with the surface disinfectant (for example Meliseptol® rapid). The concentration and the time the disinfectant must be applied, depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the device.

CAUTION!

Cleaning the Device

Do not sterilize the device.

Manufacturer's specifications

12.3 Maintenance Intervals

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. This inspection must be carried out once a year. The tests are described in Chapter Annual Inspection [► 82].

Regular inspections will assist in early detection of possible malfunctions. This helps to preserve the device and increases its safety and service life.

Two-year maintenance interval

12.4 Maintenance by Authorized Service Technician

An authorized service technician has to inspect and service the device at appropriate intervals to ensure its safety and functionality. The service interval is two years. If the service interval is not maintained, the manufacturer does not assume any liability for the functional safety of the device.

A sticker located on the rear panel of the device states the latest date for the next service or maintenance check.

Authorized service technicians are only trained and certified by the manufacturer.

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians.

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

Ask the service technician for a technical report after the service technician has inspected the device or performed any service tasks. This technical report lists the type, the scope and the results of the service as well as the date and name of the servicing company together with the signature of the service technician.

Authorized trained personnel

Unauthorized personnel

Liability

Technical documents

Technical report

WARNING!

Risk of electrical shock

To prevent electrical shock, do not open this device. Never open this device yourself. Notify the authorized service technicians of any required repairs.

WARNING!

Modification of the device

This device may not be modified without the permission of the manufacturer.

WARNING!

Modified device

If the device is modified, suitable examinations and tests must be carried out to ensure further safe use of the device.

12.5 Replacing the Fuse

CAUTION!

Replace the fuse

Before replacing the fuse, check the values for the fuse to be inserted in accordance with Chapter Technical Data [► 91].

The fuse may be defective and is in need of replacement if:

- displays and LEDs (if available on your equipment) do not light up,
- the device does not function.

Check to make sure

- the main power supply cable is properly connected to the power supply input and to a safety socket,
- the house power supply fuse is functioning.

WARNING!


Checking the fuse

Unplug the power cable from the device before checking the fuse.

Fig. 12–1 Opening the fuse holder

The device does **not** have to be opened to replace the fuse.

1. Switch device off.
2. Disconnect device from power supply by pulling the mains plug from the mains socket.
3. Remove power connection cable from mains socket.
4. The fuse holder is located next to the mains socket.
5. Remove fuse holder as depicted in Fig. Opening the fuse holder [► 76].
6. **A** Undo the latch of the fuse holder with two fingers.
7. **B** Remove the fuse holder.
8. **C** Check fuse.
9. Insert a new fuse. Use only the specified type of fuse (see Chapter Technical Data [► 91]).
10. Insert the fuse holder until it can be heard snapping into place.


12.6 Care and Maintenance of Universal Gas Connectors

The combined gas connector must be checked during the annual test and maintained at least every two years. During the maintenance, the sinter filter must be changed.

The sinter filter may be blocked or soiled due to unclean gas or gas supply residues. Heavy soiling is especially noticeable during higher flow rates. As a result, the actual gas flow is falling short of the nominal gas flow during active insufflation or the device is showing gas supply warnings although gas supply is sufficient. In this case, exchange the sinter filter.

12.6.1 Disassembly

1. Remove the connector ① using two open-end wrenches SW17, SW19.
2. Unscrew the sinter filter ② using a screwdriver if necessary.

- ① Connector
② Sinter filter

12.6.2 Accessories

Order number	Designation
Z5075-01	Connector DISS for universal gas connection (house gas supply)
Z5072-01	Sinter filter for universal gas connection
Z5077-01	Universal gas connector CO2 angled - standard

12.7 Care and Maintenance of the High-Pressure Hose

To test the high-pressure hose for leaks, do the following:

- Visual check
- Functional test

12.7.1 Visual Check

Check that the following parts are clean and not visibly damaged:

- High-pressure hose
- Connectors
- Sealing rings

12.7.2 Functional Test

The pressure drop should not exceed 5 bar after 30 minutes.

1. Connect the high-pressure hose to a standard high-pressure manometer.
2. Connect the high-pressure hose to a gas bottle.
3. Open the gas bottle and create pressure of minimum 40 bar.
4. Close the gas bottle valve.
5. Check the pressure drop after 30 minutes.

12.8 Care and Maintenance of Reusable Tube Sets (NOT FOR SALE IN USA)

Please comply with the following information when using a reusable tube set.

The reusable tube systems must be cleaned, disinfected and sterilized before each use. This applies in particular also to the first use after delivery, because all reusable tube systems are supplied non-sterile (cleaning and disinfecting after removing transport protection packing, sterilization after removing packaging). Effective cleaning and disinfection are essential prerequisites for effective sterilization. As part of your responsibility for the sterility of products during use, please make sure that

General rules

- only cleaning/disinfection and sterilization methods sufficiently validated for the device and specifically designed for the product are used.
- the used devices (RDG, sterilizer) are regularly maintained and checked.
- adherence to the validated parameters during each cycle is ensured.

Please also comply with the legal or statutory requirements applicable to your country as well as the hygiene rules of the physician's office or hospital. This applies in particular to the differing requirements regarding effective prior inactivation.

Recycling limitations

WARNING!

Cleaning of reusable tube set

The reusable tube set (NOT FOR SALE IN USA) is made of silicone, polysulfone (PSU), and stainless steel. For cleaning and disinfection, use only pH-neutral or slightly alkaline cleaning agents (e.g. neodisher MediClean 2.0 %), disinfectants (e.g. Lysetol V 8 %), and drying and rinsing agents approved for the tube set materials.

Using unsuitable agents (e.g. neodisher MediKlar rinsing agent) may damage the tube system and especially the PSU connectors.

The manufacturer has verified the reusable tube system for a limited number of recycling cycles. Always comply with the instructions on the label.

Never exceed the number of uses indicated by the manufacturer.

WARNING!

Signs of damage

Use the tear-off tabs attached to the tube to keep track of the number of sterilization cycles. The heating tube may not be sterilized after the last tab has been removed.

Please check the reusable tube (NOT FOR SALE IN USA) for signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation.

12.9 Cleaning and Disinfecting

Basics

A mechanical process (cleaning and disinfecting equipment) should be used if possible. A manual process - even when using an ultrasonic bath - should be used only if a mechanical method is not available or not possible due to the significantly reduced efficiency and reproducibility. Pretreatment is required in both cases.

Pretreatment

Process

- Rinse the products at least 1 min under running water (temperature < 35 °C/ 95 °F). Rinse all lumina of the products three times using a disposable syringe (minimum volume 100 ml).
- Place the products in a freshly prepared cleaning solution¹. With regard to selection, concentration, temperature, and exposure time, observe relevant requirements as listed in Manual cleaning and disinfecting. Use only freshly prepared solutions. Rinse all lumina of the products three times using a disposable syringe (minimum volume 100 ml).
- Manually remove all externally visible particles with a clean soft cloth and soft brush, which you are using only for this purpose.
- Rinse the products at least 1 min under running water (temperature < 35 °C/ 95 °F). Rinse all lumina of the products three times using a disposable syringe (minimum volume 100 ml).

¹If you are using a combination cleaner and disinfectant for this step - e.g., for reasons of occupational safety and health - please keep in mind that this product should be free of aldehyde (otherwise blood particles may become permanently deposited), should have a certified efficacy (e.g., VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), is suitable for the disinfection of products, and is compatible with the products being cleaned and disinfected (see Material resistance

[► 81]). Please keep in mind that the disinfectant used for the pre-treatment serves only personal protection and cannot replace the subsequent disinfection carried out after the cleaning step.

When selecting a cleaning and disinfecting equipment, make sure

- this device has a certified efficacy (e.g., DGHM or FDA approval/clearance/registration or CE marking according to DIN EN ISO 15883).
- to use a certified program for thermal disinfection (A0 value > 3000 or - for older devices - at least 5 min at 90 °C/194 °F), if possible (with chemical disinfection risk of disinfectant residues on the products).
- the used program for the products is suitable and contains sufficient rinsing cycles.
- rinsing is done only with water (e.g., purified water/highly purified water) that is sterile or low-germ (max. 10 germs/ml) and low in endotoxins (max. 0.25 endotoxin units/ml).
- the air used for drying is filtered (oil-free, low in particles and germs).
- the cleaning and disinfecting equipment will be regularly maintained and checked.

When choosing the cleaning agent system, make sure

- it is generally suitable for the cleaning of instruments made of metal and plastics.
- an additional, suitable disinfectant with certified efficacy (e.g., VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that this disinfectant is compatible with the utilized cleaning agent - unless a thermal disinfection process is applied.
- the chemicals used are compatible with the products (see Material resistance [► 81]).

It is imperative to comply with the concentrations, temperatures, and exposure times, as well as the post-rinsing requirements specified by the manufacturer of the cleaning agent and/or disinfectant.

1. Place the products into the cleaning and disinfecting equipment. Make sure the products do not touch one another. Close all lumina of the products using the existing Luer lock adapter at the rinsing port of the cleaning and disinfecting equipment.
2. Start the program.
3. Remove the products from the cleaning and disinfecting equipment after the end of the program.
4. Check and pack the products as soon as possible after removal (see Inspection [► 80], Maintenance [► 80] and Packaging [► 80] possibly after additional drying in a clean location).

Process

Documentation of the fundamental suitability of the reusable heating tube system for an effective mechanical cleaning and disinfection was prepared by an independent accredited and recognized (§15 (5) MPG) testing laboratory using the cleaning and disinfecting equipment type G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the pre-cleaning/cleaning agent neodisher MediZym (Dr. Weigert GmbH & Co. KG, Hamburg). This was done by taking into account the procedure described above.

When choosing the cleaning and disinfecting agent to be used, make sure

- they are generally suitable for the cleaning or disinfection of tube systems made of metal and plastics.
- a disinfectant with certified efficacy (e.g., VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that this disinfectant is compatible with the utilized cleaning agent.
- the chemicals used are compatible with the products (see Material resistance [► 81]).

Combined cleaning/disinfecting agents should not be used if possible. Combined cleaning/disinfecting agents can be used only in cases of very low contamination (no visible impurities).

Mechanical cleaning/disinfecting (cleaning and disinfecting equipment)

Manual cleaning and disinfecting

Cleaning procedure

It is imperative to comply with the concentrations, temperatures, and exposure times, as well as the post-rinsing requirements specified by the manufacturer of the cleaning agent and disinfectant. Use only freshly prepared solution, only water (e.g., purified water/highly purified water) that is sterile or low-germ (max. 10 germs/ml) and low in endotoxins (max. 0.25 endotoxin units/ml) and only a soft, clean and lint-free cloth and/or filtered air for drying.

1. Insert the products into the cleaning solution for the prescribed contact time and make sure the products are sufficiently submerged. Make sure the products do not touch one another.
2. Rinse all lumina of the products at least five times at the beginning or at the end of the contact time using a disposable syringe (minimum volume 100 ml). Clean the exterior with a clean, soft cloth and a soft brush that you use just for this purpose.
3. Remove the products from the cleaning solution and rinse at least three times thoroughly (for at least 1 min) with water. Rinse all lumina of the products at least five times using a disposable syringe (minimum volume 100 ml).
4. Check the products (see Inspection [► 80] and Maintenance [► 80]).

Disinfecting procedure

1. Insert the cleaned and checked products into the disinfecting solution for the prescribed contact time and make sure the products are sufficiently submerged. Make sure the products do not touch one another. Rinse all lumina of the products at least five times at the beginning or at the end of the contact time using a disposable syringe (minimum volume 100 ml).
2. Remove the products from the disinfecting solution and rinse at least five times thoroughly (at least 1 min) with water. Rinse all lumina of the products at least five times using a disposable syringe (minimum volume 100 ml).
3. Dry the products by blowing off/blowing out with filtered compressed air.
4. Package the products as soon as possible after removal (see Packaging [► 80], possibly after additional drying in a clean location).

Documentation of the fundamental suitability of the reusable heating tube system for an effective manual cleaning and disinfection was prepared by an independent accredited and recognized (§15 (5) MPG) testing laboratory using the pre-cleaning/cleaning agent neodisher Medizym and the disinfecting agent neodisher Septo 3000 (Dr. Weigert GmbH & Co. KG, Hamburg). This was done by taking into account the procedure described above.

Inspection

Check the reusable tube system after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chipping, soiling and discoloration and discard damaged products. (See label of the product for the numerical restriction of reuse). In case of damage (scratches, flaking) on the chrome-plated Luer-lock adapter, the products must be replaced immediately (no reuse permitted). Products still contaminated or dirty must be once again cleaned and disinfected.

Maintenance

Instrument oils or grease may not be used.

Packaging

Please pack the tube systems into sterile disposable packaging (single or double packaging) that comply with the following requirements (material/process):

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (temperature resistance up to min. 142 °C (288 °F) with sufficient steam permeability)
- Sufficient protection of the products or sterilization packaging against mechanical damage.

12.10 Sterilization of Reusable Tube Sets (NOT FOR SALE IN USA)

Sterilization must be carried out with one of the sterilization methods listed below; any other sterilization procedures are not permitted.

Steam sterilization

- Fractionated vacuum method ^{2,3}(with sufficient product drying⁴)
- Steam sterilizer according to DIN EN 13060/DIN EN 285 or ANSI AAMI ST79
- validated according to DIN EN ISO 17665 (valid IQ/OQ (picking) and product specific performance assessment (PQ))
- maximum sterilization temperature 138 °C (280 °F; plus tolerance according to DIN EN ISO 17665)

- Sterilization time (exposure time at sterilization temperature) at least 5 minutes at 132 °C (270 °F) / 134 °C (273 °F)

²at least three vacuum stages

³Use of the less effective gravitation procedure is allowed only in case of non-availability of the fractionated vacuum procedure; requires significantly longer sterilization times and must be validated specifically for the products, devices, and parameters in sole responsibility of the user.

⁴The actual required drying time depends on parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer condition, etc.) and must therefore be determined by the user. However, drying times should never be less than 10 minutes.

Flash sterilization procedures are not permitted. Furthermore, do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

Documentation of the general suitability of the reusable heating tube system for effective sterilization was provided by an independent accredited and recognized (§ 15 (5) MPG) testing laboratory using the autoclave Systec HX-320.

This was done by taking into account the procedure described above.

After sterilization, the products in the sterilization packaging must be stored dry and dust-free.

Storage

When selecting the cleaning and disinfecting agents, please ensure that the following ingredients are not included:

Material resistance

- organic, mineral and oxidizing acids (minimum admissible pH 6.5)
- strong alkalis (maximum admissible pH 9.5, enzymatic/neutral or slightly alkaline cleaner recommended)
- organic solvents (alcohols, ethers, ketones, petroleum)
- oxidizing agents (such as hydrogen peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

Never clean the reusable tube with metal brushes or steel wool. Rinse aids or acid neutralizers are not permitted. All products may not be exposed to temperatures in excess of 142 °C (288 °F)! Failure to comply with these instructions voids any liability.

Measured values and tolerances

13 Annual Inspection

Each test conducted has to be documented with date and signature on the test log. The following measuring tools and resources were used by the manufacturer to determine the listed measurements and tolerances:

Manometer	Range 0-100 mmHg, error class 1.6
Syringe	60 ml
Silicone tube	8 mm x 2 m
T adapter	8-8-8 mm
Veress cannula	length 100 mm opening diameter 1.4 mm, inner cannula diameter 1.6 mm

An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

13.1 Electrical Safety Test

- Perform a visual inspection. Make sure that
 - the fuse corresponds with the specifications indicated by the manufacturer,
 - labels and stickers on device are legible,
 - the mechanical condition of the device allows for its safe use,
 - the device is clean to ensure proper and safe functionality.
- Carry out the measurements for the earth leakage current, short-circuit current/equipment leakage current and protective earth resistance by using an accessible conductive part as per IEC 62353 in the current version or according to the applicable national standard.
- If specific electrical safety tests are required for this device, please refer to the corresponding section below in this chapter. Only required specific electrical safety tests are described.

13.2 Basic Function Test (in Standard Operating Mode)

- Remove insufflation tube from device.
- Use the ON/OFF button (push-push) to turn the device on. The device now conducts a device check. A short beep can be heard. Set the Standard operating mode.
- The factory default settings are 15 mmHg for the nominal pressure and 3 l/min for the nominal flow.
- The display depicts the following values:
 - Nominal pressure 15* [mmHg]
 - Nominal gas flow 3* [l/min]
 - Actual pressure 0 [mmHg]
 - Gas consumption 0.0 [l]

*These values correspond with the factory setting. If values in the configuration menu were changed, these changed values are displayed.
- Start insufflation: Press the function field **Start**. The display depicts the following values:
 - Actual pressure 0 [mmHg]

Veress insufflation

Streaming gas can be heard at the insufflation tube.
- Select the maximal nominal gas flow. The display depicts the following values:
 - Nominal gas flow maximal value [l/min]
 - Actual pressure 0 [mmHg]

Insufflation: Standard

Streaming gas can be heard at the insufflation tube.

7. Stop insufflation: Press the function field **Stop**.

The display depicts the following values:

Actual pressure 0.0 [mmHg]


Gas consumption > 0.0 [l]

8. Press the gas volume function field.

Gas consumption 0.0 [l]

The basic function check of the device is complete.

13.3 Testing the Pressure Sensors in Standard Mode Operating

1. Select **Standard operating mode**.
2. Select a nominal gas flow rate of 1.0 l/min.
Do not press the **Start/Stop** function field.
3. Connect a manometer (1) and an air-filled syringe (2) to the insufflation tube connection (3)/adapter (4).
4. Use the syringe to generate a pressure of at least 10 mmHg, which is indicated on the manometer.
Actual pressure display: 10 ± 2 [mmHg]
5. Use the syringe to generate a pressure of at least 20 mmHg, which is indicated on the manometer.
Actual pressure display: 20 ± 2 [mmHg]
6. Use the syringe to generate a pressure of at least 30 mmHg, which is indicated on the manometer.
Actual pressure display: 30 ± 2 [mmHg]

CAUTION!

Gas extraction

Never use the syringe to extract gas from the device.

13.4 Pressure Monitoring Test (in Standard Operation Mode)


1. Set the **Standard operating mode**.
2. Select a nominal pressure of 15 mmHg and a nominal gas flow of 3 l/min.
3. In the configuration menu **Venting valve pressure** threshold, set a venting pressure of 3 mmHg.
4. Use the syringe to generate a pressure of at least 19 mmHg, which is indicated on the manometer. Start insufflation: Press the function field **Start**. An acoustic warning signal is emitted with a pressure of more than 19 mmHg (for 5 seconds) and the display depicts **Overpressure**.
5. Reduce the pressure. The warning ends when the pressure falls below 19 mmHg (nominal pressure plus 4 mmHg). Stop insufflation: Press the function field **Stop**.
6. Select a nominal pressure of 29 mmHg.
7. Use the syringe to generate a pressure of at least 30 mmHg, which is indicated on the manometer. Start insufflation: Press the function field **Start**. An acoustic warning sound (beep) is emitted without delay in case pressure exceeds 30 mmHg and the display depicts **Overpressure**.
8. Reduce the pressure. The warning ends when the pressure falls below 30 mmHg. Stop insufflation: Press the function field **Stop**.

13.5 Venting Valve Test

See also Chapter Testing the Pressure Sensors in Standard Mode Operating [► 83].

1. In the configuration menu **Venting system** set the venting system to ON (with setting of **Veress insufflation OFF**).
2. In the configuration menu **Venting valve response time** set a venting time of 3 seconds.
3. In the configuration menu **Venting valve pressure threshold**, set a venting pressure of 3 mmHg.
4. Select a nominal pressure of 15 mmHg and a nominal gas flow of 10 l/min.
5. Use the syringe to generate a pressure of at least 18 mmHg, which is indicated on the manometer. Start insufflation. With a pressure of more than 18 mmHg (for 3 seconds), the venting valve is activated and the display depicts **Venting system active**.

13.6 Maximal Device Pressure Test

1. Select **Standard operating mode**.
2. Select the maximal nominal gas flow.
3. Connect a manometer (1) and an open Veress cannula (2) to the insufflation tube connection (3)/adapter (4).
4. Select the maximal gas flow.
5. Start insufflation:
Press the **Start** function field. The manometer registers a pulsing pressure increase. When the pressure stabilizes, the manometer registers a maximum pressure of between 55 and 65 mmHg (Bariatric mode: 65 to 75 mmHg).
6. Stop insufflation:
Press the **Stop** function field.

13.7 Gas Flow Rate Test

Test setup with open connection, without connected insufflation tube.

- Select a nominal gas flow rate of 15 l/min.
 - Start insufflation:
Press the **Start** function field.
 - Press the gas volume function field (0.0 l must be displayed).
 - Now start measuring for one minute.
 - Stop the insufflation after one minute: Press the **Stop** function field.
- The gas consumption should be at least 11-12 l.

Each successfully conducted test must be documented with the test log.

14 Electromagnetic compatibility

CAUTION!

Accessories

To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device **Insufflator 50L** must be used only with the accessories listed in Chapter Accessories [► 92].

Precautionary measures

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

The device is to be used only for the purposes described in the Instructions for Use and is intended for use in environments in Professional Healthcare Facility Environment. This applies even if individual requirements meet the conditions for deviating electromagnetic environments. During installation and commissioning as well as during operation of the device, the compliance with the notes and instructions for EMC must be strictly observed.

To ensure the basic safety and essential functionality in relation to electromagnetic interference over the life of the device, the device must be restarted after 24 hours so that a diagnostic self-test can be performed. The maintenance intervals indicated in Chapter Maintenance Intervals [► 74] must also be observed.

This device complies with the electromagnetic compatibility (EMC) requirements for medical electrical devices as defined by IEC 60601-1-2. The limits used in testing provide a basic level of safety against typical electromagnetic interference likely to occur in professional health care facilities. Nevertheless, it can happen that individual performance features are no longer available or only to a limited extent due to the presence of EM interference.

In the presence of EM interference, it may happen that the heating function is affected. If the heating function is influenced by EM interference, it will be switched off for safety reasons. After the EM interference has faded, the heating function is reactivated.

Furthermore, it is possible that interruptions of the supply voltages causes the insufflation being stopped because safe operation is not guaranteed when falling below the specified input voltage. After restoration of the specified input voltage, insufflation is automatically restarted.

In some cases of long voltage interruptions the device can be reset. Pressure in cavity is held in that time.

14.1 Electrical Connections

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

Hospital employees should be informed about and trained in ESD precautionary measures.

14.2 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The device **Insufflator 50L** is intended for use in an environment as described below. The user/operator of the device **Insufflator 50L** should make sure the device is operated within such an environment.

ESD (Electrostatic Discharge) precautionary measures

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device Insufflator 50L uses RF energy solely for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device Insufflator 50L is suitable for use in all establishments, other than public establishments and those directly connected to the public mains network that supplies buildings used for public purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

14.3 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity

The device **Insufflator 50L** is intended for use in an electromagnetic environment as described below. The user/operator of the device **Insufflator 50L** should make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	IEC 60601 test level	Compliance levels	Electromagnetic environment guidelines
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV contact ± 2, ± 4, ± 8, ± 15 kV air	± 8 kV contact ± 2, ± 4, ± 8, ± 15 kV air	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/bursts according to IEC 61000-4-4	± 2 kV for power lines Repetition frequency 100 kHz	± 2 kV for power lines Repetition frequency 100 kHz	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Surges according to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brownouts, and fluctuations of the power supply according to IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	
	0 % UT; 250/300 cycles	0 % UT; 250/300 cycles	
Supply frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.
Conducted RF IEC 61000-4-6	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Transients RF IEC 61000-4-3	3 V/m 0.15 MHz to 80 MHz 80 % AM by 1 kHz	3 V/m 0.15 MHz to 80 MHz 80 % AM by 1 kHz	

IMMUNITY to proximity fields from RF wireless communications equipment**Test method IEC 61000-4-3**

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sinus	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

**WARNING!****Portable HF communication equipment**

Portable HF communication equipment can affect the performance characteristics of the device Insufflator 50L. Such equipment must therefore comply with a minimum distance of 30 cm (regardless of all calculations) from the device Insufflator 50L, its components and cables.

15 Error and Warning Messages

Error and warning messages	Cause	Troubleshooting
Check gas supply	(During device check) Existing gas supply pressure is too low.	• Open gas bottle or • Replace the gas bottle. • Check the house gas supply.
	(During surgery) The gas supply pressure has dropped below 15 bar.	• Detach the insufflation tube. • Close gas supply valve. • Replace gas supply. • Open gas supply valve. • Connect insufflation tube.
	Insufficient house gas supply.	• Check the house gas supply. • Open or remove pressure reduction valve if in line with gas bottle. • Check if appropriate house gas connector is used.
Change gas bottle!	The gas supply pressure has dropped below 30 bar.	• Prepare for changing the gas bottle.
Overpressure	The pressure monitor shows that the actual pressure is at least 2-5 mmHg (depending on the settings in the configuration menu) above the nominal pressure.	• Determine the cause for exceeding the nominal pressure. Check the electronic controls of the device if overpressure exists for a longer period of time (see Chapter Annual Inspection [► 82]).
	The actual pressure has reached 30 mmHg/ 20 mmHg (depending on the mode).	• Reduce the nominal pressure and determine the cause of exceeding the nominal pressure if necessary.
Venting system active	The actual pressure is at least 2-5 mmHg and 2-5 seconds above the nominal pressure.	• Determine the cause for exceeding the nominal pressure. Check the electronic controls of the device if overpressure exists for a longer period of time (see Chapter Annual Inspection [► 82]).
Overpressure/venting system active.	The pressure monitor shows that the actual pressure is 2-5 mmHg for 2-5 seconds above the nominal pressure. The overpressure was not reduced within 5 seconds by the venting system.	• Determine the cause for exceeding the nominal pressure. Check the electronic controls of the device if overpressure exists for a longer period of time (see Chapter Annual Inspection [► 82]). Reduce the nominal pressure. Check if the instrument's stopcock is open or the tube is plugged up.
Occlusion	Tube or instrument occlusion.	• Localize the cause and open/eliminate the occlusion.
	Faulty Veress needle insertion. Stopcock is closed.	• Check that the Veress cannula is positioned correctly in the abdomen and make sure the instrument's stopcock is open.
Contamination/Call service	The device is contaminated with fluid.	• The message is repeated with each Start/Stop. It is possible to continue using the device with this error message until the device is turned off with the ON/OFF button (push-push). • The device must be checked by an authorized service technician or clearly marked with a label referring to the contamination and then twice enclosed in a safety foil, sealed, and returned to the manufacturer for repairs.

Error and warning messages	Cause	Troubleshooting
Gas heater defective/Call service	Gas heater malfunction.	Check gas heating with a new tube. If this error message is displayed again, have the gas heating checked by a qualified service technician. It is possible to continue using the device without gas heating.
Gas temperature >42 °C The device terminates gas heating	The temperature of the gas exceeds 42 °C.	- The device insufflates without heating the gas. - Continue the operation without gas heating or use a different gas heating tube. - Turn device off and on again. This reactivates the gas heating. - Check gas heater after surgery using a different tube. If this error message is displayed again, have the gas heating checked by a qualified service technician. It is possible to continue using the device without gas heating.
Error message/Call service: Electronic defective Sensor defective Calibration error	The device does not work properly The internal safety system has been triggered.	- Switch the device off and back on after approx. 3 seconds have expired using the ON/OFF button (push-push). - The device is defective if the error message is displayed again. Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.
Device temperature error/Turn device off Insufflation is stopped or cannot be started.	The device temperature is above 70 °C.	Use the ON/OFF button (push-push) to turn the device off and allow it to cool off for about 10 minutes. Make sure the device is not set up in the vicinity of heat sources.
	The device temperature is below 10 °C.	Turn the device off for about 10 minutes using the ON/OFF button (push-push). The surrounding room temperature must be above 10 °C.
Venting valve defective	Venting system malfunctioning.	Use the ON/OFF button (push-push) to turn the device off and turn device back on after approx. 3 seconds. If this error message is displayed again after the device check has concluded, have the device checked by an authorized service technician. The device can be still operated but without venting system. The error message is repeated with each Start/Stop.

16 Technical Data

Power connection:	100-240 V~ Mains fuse T 3.15 A H 250 V Connection for potential equalization
Frequency:	50/60 Hz
Maximal power consumption:	150 VA
Current consumption:	100 V~ 60 Hz: 0.45 A (US, UK, CA: 0.45 A) 240 V~ 50 Hz: 0.25 A (US, UK, CA: 0.25 A)
Protection class:	I
Protection level:	Type BF
Moisture protection:	IP21
Dimensions:	Width x Height x Depth 380 x 142 x 432 (mm)
Weight:	approx. 10 kg/22 lbs
Operating conditions:	10 to 40 °C / 50 to 104 °F 30 to 75 % rel. humidity 70 to 106 kPa ambient pressure maximum usage altitude: 3000 m above sea level
Transportation conditions:	-20 to +60 °C / -4 to +140 °F 10 to 85 % rel. humidity 70 to 106 kPa air pressure
Storage conditions:	-20 to +60 °C / -4 to +140 °F 10 to 85 % rel. humidity 70 to 106 kPa air pressure
Insufflation medium:	Medical grade CO ₂
Maximum output pressure:	75 mmHg (1 mmHg = 1.33 mbar = 133 Pa)
Maximum gas supply pressure:	80 bar/1160 PSI
Minimum gas supply pressure (bottle):	15 bar/218 PSI
Minimum gas supply pressure (house gas):	3.4 bar/50 PSI
Maximum gas flow:	Bariatric: 50 l/min Standard: 40 l/min Pediatric: 20 l/min Vessel Harvest: 10 l/min
Pressure range:	1-20 mmHg (Pediatric/Vessel Harvest) 1-30 mmHg (Standard/Bariatric)
Pressure display resolution:	1 mmHg
Gas flow display resolution:	1 l/min (Standard/Bariatric/Vessel Harvest/Pediatric in range of 2 to 20 l/min) 0.1 l/min (Pediatric in range of 0 to 1.9 l/min)
Volume display resolution:	in range of 0 to 9.9 l resolution = 0.1 l in range of 10 to 999 l resolution = 1 l
Service connections/ports:	Service interface
Hose for central gas supply:	Design version: Connections are stainless steel and chrome-plated brass, PVC tubing, gray Output Permissible operating pressure: 1400 kPa (14 bar) Burst pressure: > 5600 kPa (56 bar) at 23 °C Max. flow for compressed gases > 100 l/min at 320 kPa Wall connection: Gas type specific plug inserts according to different national standards Device-side connection: SAE J512, 45° Max operating pressure 80 bar
High pressure hose:	

17 Accessories

Order number	Designation
Disposable Insufflation Tube Sets	
Z0293	Tube set for insufflation, disposable
Z0260	Heating tube set for insufflation, disposable
Reusable Insufflation Tube Sets and Accessories (NOT FOR SALE IN THE USA)	
Z0452-01	Tube set for insufflation, reusable (20 times)
Z0280	Heating tube set for insufflation, reusable (100 times)
Z0536-01	Filter, single-use
Z0164-01	ISO adapter
Optional Accessories	
Z0531-00	Switch valve for insufflators
Z5075-01	DISS adapter for universal gas connector
Z5077-01	Universal gas connector CO ₂ angled - standard
Hoses for Gas Supply	
Z5044-01	Hose for CO ₂ gas bottle DIN:US, 1.5 m
Z5045-01	Hose for CO ₂ gas bottle ISO:US, 1.5 m
Z0499-01	Hose for CO ₂ gas bottle Pin:US, 1.5 m
Z5080-01	Hose for CO ₂ gas bottle CGA:US, 1.5 m
Z5100-01	Hose for CO ₂ central gas DIN:US, 3 m
Z5101-01	Hose for CO ₂ central gas DIN:US, 5 m
Z5102-01	Hose for CO ₂ central gas NF:US, 3 m
Z5103-01	Hose for CO ₂ central gas NF:US, 5 m
Z5104-01	Hose for CO ₂ central gas UNI:US, 3 m
Z5105-01	Hose for CO ₂ central gas UNI:US, 5 m
Z5106-01	Hose for CO ₂ central gas AGA:US, 3 m
Z5107-01	Hose for CO ₂ central gas AGA:US, 5 m
Z5108-01	Hose for CO ₂ central gas DISS:DISS, 3 m
Miscellaneous	
Z0102-01	Power cord US, 2.5 m
Z0101-01	Power cord Euro, 2 m
Z0176-00	Power cord UK, 2 m
Z0178-01	Power cord China, 1.8 m
Z0560-01	Power cord Brasil, 1.8 m
Z5040-01	Fuses T3.15A, 5x20 mm (10 pieces)
Z5072-01	Sinter filter for universal gas connector (5 pieces)

Please fill out this form when returning the device:

Name of owner:

--

Sales partner:

Address of person returning unit:

Street:

House number:

--	--

ZIP/Postal code:

City:

--	--

Country:

IMPORTANT!

Serial number (see identification plate):

Device type:

Description of defect:

[illegible]

Contact

Signature

Date _____

Glossary

EMC

The electromagnetic compatibility describes the ability of a device to function satisfactorily within an electromagnetic environment without adding unacceptable electromagnetic interferences/disturbances to the environment that may cause problems for other devices or equipment located nearby.

ME Device

Medical electrical device used for therapy, monitoring or diagnosis of patients, equipped with no more than one connection to a supply network and which necessarily comes into physical or electrical contact with the patient or which transfers energy to or from the patient or which records or captures such energy transfer to or from the patient.

ME System

Medical electrical system consisting of a combination of devices, of which at least one is classified as a medical electrical device and specified by the manufacturer as such and which are joined together by a functional connection or by using a power strip.

Index

A

Authorized service technician	10
Authorized trained personnel	75
Automatic venting system	72

B

Basics	78
--------	----

C

Care and maintenance	10
Characteristics	36, 43, 51, 59
Cleaning procedure	80
Configuration menu II	69
Contamination	10
Contamination/Call service	72
Contraindications	43, 59

D

Device Check	39, 47, 55
Disconnect device from power supply	28, 35
Disinfecting procedure	80

E

Error message Contamination	72
Error message when starting with defective venting valve	73
Essential Performance	11
Exceeding nominal pressure	72
Exceeding nominal pressure for more than 5 seconds	72
Exchange of the gas supply adapter	28
Exchangeable sinter filter	26
Exclusion of liability	10

F

Federal Law (only for U.S. market)	10
Flow safety limit	45

G

Gas supply displays	30
Gas supply with gas bottle	30
Grounding contact	25

H

High Pressure Hose Types	27
House gas supply	31

I

Increasing/decreasing nominal pressure	38
Inspection	80
Installation	27
Insufflation parameters	36, 43, 51, 59
Insufflation source	72
Insufflation: Standard	38

L

Liability	75
-----------	----

M

Mains connection	25
Maintenance	80
Manual cleaning and disinfecting	79
Manufacturer's specifications	74
Material resistance	81
Mechanical cleaning/disinfecting	79

N

Nominal pressure limit 30 mmHg	72
--------------------------------	----

O

Only for U.S. operators	25
-------------------------	----

P

Packaging	80
Patients	36, 43, 51, 59
Potential equalization	25
Power failure	73
Preparation	39, 46, 61
Pressure safety limit	38
Pretreatment	78
Process	79

R

Reset the display	35
Returning the device	24

S

Setting the venting system response time	68
Setting the venting valve pressure limit	67
Setting the venting valve status	67
Setting up	24
Steam sterilization	80
Storage	81
Subject to technical changes	9

T

Technical documents	75
Technical report	75
Testing gas heating functionality	40, 56
Text message - >Call service	73
Two-year maintenance interval	74

U

Unauthorized personnel	75
Universal gas connector	25

V

Veress insufflation	38
---------------------	----

W

Warning message Occlusion	72
Waste management	10
With deactivated venting system or during Veress insufflation	72

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